



Disposable Bed Pads



Fixing and Maintenance Instructions

INTRODUCTION

Thank you for deciding to purchase the Disposable Bed Pads from Aidapt. These products are manufactured from the finest material for reliable, trouble-free service.

NB. This equipment must be installed by a competent person.

Code	Size (mm)	SAP	No. in Pack 25
VM845A	600x600	3	
VM845AB	600x600	5	25
VM848A	900x600	7	25

BEFORE USE

Carefully remove all packaging. Avoid using any knives or other sharp instruments as this may damage the surface of the product.

Check the product for any visible damage. If you see any damage or suspect a fault, do not use your product, but contact your supplier for support.

INTENDED USE

Aidapt's Disposable Bed Pads are disposable incontinence pads for men and women, designed to offer a hygienic and comfortable solution to protect beds and chairs.

With a soft facing material for comfort and a waterproof backing, Aidapt's incontinence pads are available in different sizes and absorbency rates.

Absorbency rates for incontinence pads are given as an SAP (Super Absorbent Polymers) rating; a higher rating shows a greater level of absorbency.

FITTING YOUR DISPOSABLE BED PADS Simply

place the Disposable Bed Pad onto the chair, bed or other surface, with the waterproof backing facing down. No fasteners or adhesives are required.

NB. You can place the Disposable Bed Pad directly onto the mattress and fit a sheet over it, or position it above the bed sheet, as preferred.

DISPOSAL

This product is disposable and no attempt should be made to clean and re-use it. Change the Disposable Bed Pad as often as required, this will vary depending on the user and environment. Dispose of the used Bed Pad in a plastic bag with your non-recyclable rubbish.

CARE & MAINTENANCE

Please make a safety check of the product at regular intervals or if you have any concerns.

IMPORTANT INFORMATION The

information given in this instruction booklet must not be taken as forming part of or establishing any contractual or other commitment by Aidapt Bathrooms Limited, Aidapt (Wales) Ltd or its agents or its subsidiaries and no warranty or representation concerning the information is given.

Please exercise common sense and do not take any unnecessary risks when using this product; as the user you must accept liability for safety when using the product.

Please do not hesitate to contact the person who issued this product to you or the manufacturer (detailed below) if you have any queries concerning the assembly/use of your product.

Aidapt Bathrooms Ltd, Lancots Lane, Sutton Oak, St Helens, WA9 3EX Telephone: +44 (0) 1744 745 020 • Fax: +44 (0) 1744 745 001 • Web: www.aidapt.co.uk Email: accounts@aidapt.co.uk • adaptations@aidapt.co.uk • sales@aidapt.co.uk



Disposable mattress pads



Use and Maintenance Instructions

INTRODUCTION

Thank you for purchasing disposable mattress pads from Aidapt. These are made with the best materials available. If you use them correctly, they will be of great service to you.

NB This equipment must be installed by a competent person.

Code	Dimensions (mm)	SAP	Quantity of units in Packaging 25
VM845A	600x600	3	
VM845AB	600x600	5	25
VM848A	900x600	7	25

BEFORE USE

Check that the product is in good condition before using it. If you notice any damage or suspicious defect, do not use the product and contact your dealer. Store the packaging carefully.

RECOMMENDED USE Aidapt

disposable pads have been designed for men and women suffering from incontinence.

They are intended to provide you with a hygienic and comfortable solution and are intended to protect beds and chairs.

The front is made of soft materials for your comfort. The back side is waterproof.

Aidapt mattress pads are available in different sizes and different absorption rates.

The absorption rates for pads intended for people suffering from incontinence are expressed in SAP (Super Absorbent Polymer) rates; the higher the rate, the greater the absorption capacity.

Delivered in packs of 15 or 20.

INSTALLING DISPOSABLE PADS Simply place

the disposable pad on the chair, bed or other surface you wish to protect, with the waterproof side down. No fasteners or adhesives are required.

NB: You can place the disposable mattress pad directly on the mattress then place a sheet on top or place the mattress pad on the sheet, depending on your preference.

DISPOSAL This

product is disposable. It cannot be washed or reused. Change the pad as often as necessary, this varies depending on the user and environment. Dispose of the mattress pad in a plastic trash bag with non-recyclable trash.

INTERVIEW

Check your product regularly or in case of doubt.

IMPORTANT INFORMATION The

information given in this user guide should not be construed as forming part of or establishing any contract or commitment by Aidapt Bathrooms Limited Aidapt Wales nor its agents or subsidiaries and no warranty is given in relation to this information.

Please use common sense and not take unnecessary risks when using the product; as a user, you must accept your share of responsibility when using the product.

For further information, please contact your reseller.

Altai France (Altai Europe Ltd)

ZA du Bel Air, 17D rue de Creil, 95340 BERNES-SUR-OISE, France Tel: 01.30.34.90.84

| Fax: 01.30.34.00.68 | E-mail: altaieurope@orange.fr



Disposable pad



Assembly and maintenance instructions

INTRODUCTION

Thank you for purchasing this Aidapt disposable pad. This product is made from the best materials available.

Danger! This device must be installed by qualified personnel!

Code	Size (mm) 600x600	SAP	Number per pack 25
VM845A		3	
VM845AB	600x600	5	25
VM848A	900x600	7	25

BEFORE USE

Carefully remove the packaging. Do not use knives or other sharp instruments as they may damage the surface of the product.

Please check your product for visible damage before use. If you notice any damage or suspect a fault, please do not use the product and instead contact your dealer.

INTENDED USE Aidapt's

disposable underpads are single-use incontinence mats for men and women. They are designed to provide a hygienic and comfortable solution for protecting beds and chairs.

Aidapt's pads have a soft upper for comfort and a waterproof side. The incontinence mats are available in different sizes and with different absorption levels.

Absorption strengths for incontinence mats are measured in SAP (Super Absorbent Polymers); a larger value means a higher absorption strength.

ATTACHING THE PAD Simply place

the pad, waterproof side down, on a chair, bed, or other surface. Neither fastening aids nor adhesives are required.

Note: You can either place the pad on the mattress and put a sheet over it, or place it directly on the sheet, depending on your preference.

DISPOSAL This

product is intended for single use and no attempt should be made to clean and reuse the product. Change the pad as often as necessary; this depends on the user and the situation. Dispose of the pad in a plastic bag and put it in the non-recyclable trash.

CARE & MAINTENANCE

Please perform a safety check of the product periodically or if you have any concerns.

IMPORTANT INFORMATION The

information provided in this instruction manual does not form the basis of or form part of any contractual or other obligation of Aidapt Bathrooms Limited, Aidapt (Wales) Ltd. or any of its agents or its affiliates and no warranty or binding representation is made with respect to the information in question.

Please use common sense and do not take unnecessary risks when using this product; As a user, you are responsible for the safety when using the product.

Please do not hesitate to contact the person who issued this product to you or the manufacturer (described below) if you have any questions regarding assembly/use of your product.

Altai Germany (Altai Europe Ltd.)
Complaints & returns to: WH
Elektronik, Merowingerstraße 83, 47533 Kleve Tel:
02821/971620 | Fax: 02821/12815