



EC Declaration Of Conformity



Physio-Control, Inc.
11811 Willows Road NE
Redmond, WA 98052 USA



**Stryker European
Operations B.V.**
Herikerbergweg 110
1101 CM Amsterdam
The Netherlands

Notified Body 0123
TÜV SÜD Product Service GmbH
Ridlerstrasse 65
80339 München, Germany

PHYSIO-CONTROL declares that the CE-marked product

ITEM
LIFEPAK® CR2 Defibrillator

PART NUMBER
See Table 1

Conforms to European Community Council Directive 93/42/EEC (Medical Device Directive), as amended through 2007/47/EC, and is a Class IIb Device assessed under Annex II.

The CE marked product also conforms to the following EC Directives:

DIRECTIVE

2012/19/EU
2014/53/EU
2011/65/EU

SUBJECT

Waste Electrical and Electronic Equipment Directive
Radio Equipment Directive (RED)
Restriction of the use of certain hazardous substances (RoHS Directive)
• Annex I – Category 8 – Medical Devices
• Annex IV – Exemption 17

The CE marked product complies with the following standards:

STANDARD

EN 60601-1:2006+A12:2014
EN 60601-1-2:2015
EN 60601-1-6: 2010+A1:2015
EN 62366:2008+A1:2015
EN 60601-1-11:2015

EN 60601-2-4:2011
EN 62304:2006+A1:2015
EN1041:2008+A1:2013
EN ISO 15223-1:2012

EN 50419:2006

EN ISO 10993-1:2009+AC: 2010
EN ISO 14971:2012
EN 62281:2013

SUBJECT

General requirements for safety for medical electrical equipment
Electromagnetic disturbances – Requirements and tests
Safety requirements for usability
Application of usability engineering to medical devices
Medical electrical equipment Part 1-11: General requirements for basic safety and essential performance - Collateral standard: Requirements for medical electrical equipment and medical electrical systems used in the home healthcare environment
Safety requirements for cardiac defibrillators
Medical Device Software – Software Life Cycle Processes:
Information supplied by the manufacturer with medical devices
Symbols to be used with medical device labels, labeling and information to be supplied
Marking of electrical and electronic equipment in accordance with Article 11(2) of Directive 2002/96/EC (WEEE)
Biological evaluation of Medical Devices Part 1: Evaluation and testing
Application of risk management to medical devices
Safety of primary and secondary lithium cells and batteries during transport

The CE marked product was evaluated with the following accessories:

LPCR2 Carry Case

LPCR2 Handle

Separately CE Marked accessories:

QUIK-STEP™ Pacing/ECG/Defibrillation Electrodes
Ambu First Responder Kit

LIFEPAK CR2 Lithium Manganese Dioxide Battery Pack

Non-Medical accessories:

LPCR2, Trainer
Carry Tote, Trainer
USB2.0 Male to Micron B Cable

Table 1

Ref #	Part Number	Description	GMDN
99512-000064	CR2-2-000069	LPCR2,CELL,B,NL,DE,S,1,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,360J,360J,50J,75J,90J,ALRT,NL	48047
99512-000065	CR2-2-000070	LPCR2,WIFI,B,NO,EN-GB,A,2,ROS,MODET,ROS,FLX,CKB,CPR180,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,NO	48047
99512-000066	CR2-2-000071	LPCR2,WIFI,H,NO,EN-GB,A,2,ROS,MODET,ROS,FLX,CKB,CPR180,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,NO	48047
99512-000067	CR2-2-000072	LPCR2,CELL,B,NO,EN-GB,A,2,ROS,MODET,ROS,FLX,CKB,CPR180,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,NO	48047
99512-000068	CR2-2-000073	LPCR2,CELL,H,NO,EN-GB,A,2,ROS,MODET,ROS,FLX,CKB,CPR180,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,NO	48047
99512-000069	CR2-2-000074	LPCR2,WIFI,B,FI,EN-GB,A,2,ROS,MODET,ROS,FLX,CKB,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,FI	48047
99512-000070	CR2-2-000075	LPCR2,WIFI,H,FI,EN-GB,A,2,ROS,MODET,ROS,FLX,CKB,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,FI	48047
99512-000071	CR2-2-000076	LPCR2,CELL,B,FI,EN-GB,A,2,ROS,MODET,ROS,FLX,CKB,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,FI	48047
99512-000072	CR2-2-000077	LPCR2,CELL,H,FI,EN-GB,A,2,ROS,MODET,ROS,FLX,CKB,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,FI	48047
99512-000073	CR2-2-000078	LPCR2,WIFI,H,NL,EN-GB,S,2,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,NL	48047
99512-000074	CR2-2-000079	LPCR2,WIFI,H,NL,EN-GB,A,2,ROS,MODET,ROS,FLX,CKB,CPR180,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,NL	48047
99512-000075	CR2-2-000080	LPCR2,WIFI,B,DA,EN-GB,A,2,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,DK	48047
99512-000076	CR2-2-000081	LPCR2,WIFI,H,DA,EN-GB,A,2,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,DK	48047
99512-000077	CR2-2-000082	LPCR2,CELL,B,DA,EN-GB,A,2,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,DK	48047
99512-000078	CR2-2-000083	LPCR2,CELL,H,DA,EN-GB,A,2,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,DK	48047
99512-000079	CR2-2-000084	LPCR2,WIFI,H,NL,FR,S,2,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,BE	48047
99512-000080	CR2-2-000085	LPCR2,WIFI,H,NL,FR,A,2,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,BE	48047
99512-000081	CR2-2-000086	LPCR2,WIFI,H,FR,NL,S,2,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,BE	48047
99512-000082	CR2-2-000087	LPCR2,WIFI,H,FR,NL,A,2,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,BE	48047
99512-000083	CR2-2-000088	LPCR2,WIFI,H,FR,DE,S,2,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,LU	48047
99512-000084	CR2-2-000089	LPCR2,WIFI,H,FR,DE,A,2,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,LU	48047
99512-000085	CR2-2-000090	LPCR2,WIFI,B,SV,EN-GB,A,2,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,SE	48047
99512-000086	CR2-2-000091	LPCR2,WIFI,H,SV,EN-GB,A,2,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,SE	48047
99512-000087	CR2-2-000092	LPCR2,CELL,B,SV,EN-GB,A,2,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,SE	48047
99512-000088	CR2-2-000093	LPCR2,CELL,H,SV,EN-GB,A,2,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,SE	48047
99512-000055	CR2-2-000094	LPCR2,WIFI,B,IT,EN-GB,A,2,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,CH	48047
99512-000056	CR2-2-000095	LPCR2,WIFI,H,IT,EN-GB,A,2,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,CH	48047
99512-000057	CR2-2-000096	LPCR2,WIFI,B,FR,EN-GB,A,2,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,CH	48047
99512-000058	CR2-2-000097	LPCR2,WIFI,H,FR,EN-GB,A,2,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,CH	48047
99512-000059	CR2-2-000098	LPCR2,WIFI,B,DE,EN-GB,A,2,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,CH	48047
99512-000060	CR2-2-000099	LPCR2,WIFI,H,DE,EN-GB,A,2,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,CH	48047
99512-000061	CR2-2-000100	LPCR2,WIFI,B,EN-GB,DE,S,1,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,IE	48047
99512-000062	CR2-2-000101	LPCR2,WIFI,H,EN-GB,DE,S,1,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,IE	48047
99512-000063	CR2-2-000102	LPCR2,WIFI,B,EN-GB,DE,A,1,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,IE	48047
99512-000089	CR2-2-000103	LPCR2,WIFI,H,EN-GB,DE,A,1,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,IE	48047
99512-000090	CR2-2-000104	LPCR2,WIFI,B,EN-GB,DE,S,1,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,GB	48047
99512-000091	CR2-2-000105	LPCR2,WIFI,H,EN-GB,DE,S,1,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,GB	48047
99512-000092	CR2-2-000106	LPCR2,WIFI,B,EN-GB,DE,A,1,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,GB	48047
99512-000093	CR2-2-000107	LPCR2,WIFI,H,EN-GB,DE,A,1,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,GB	48047
99512-000094	CR2-2-000108	LPCR2,WIFI,B,DE,EN-GB,A,1,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,DE	48047
99512-000095	CR2-2-000109	LPCR2,WIFI,H,DE,EN-GB,A,1,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,DE	48047
99512-000096	CR2-2-000110	LPCR2,WIFI,B,DE,EN-GB,A,1,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,DE	48047

Ref #	Part Number	Description	GMDN
99512-000097	CR2-2-000111	LPCR2,WIFI,H,DE,EN-GB,A,1,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,DE	48047
99512-000098	CR2-2-000112	LPCR2,WIFI,B,DE,EN-GB,A,2,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,DE	48047
99512-000099	CR2-2-000113	LPCR2,WIFI,H,DE,EN-GB,A,2,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,DE	48047
99512-000100	CR2-2-000114	LPCR2,WIFI,B,DE,EN-GB,A,1,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,AT	48047
99512-000101	CR2-2-000115	LPCR2,WIFI,H,DE,EN-GB,A,1,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,AT	48047
99512-000102	CR2-2-000116	LPCR2,WIFI,B,DE,EN-GB,A,1,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,AT	48047
99512-000103	CR2-2-000117	LPCR2,WIFI,H,DE,EN-GB,A,1,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,AT	48047
99512-000104	CR2-2-000118	LPCR2,WIFI,B,DE,EN-GB,A,2,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,AT	48047
99512-000105	CR2-2-000119	LPCR2,WIFI,H,DE,EN-GB,A,2,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,AT	48047
99512-000106	CR2-2-000120	LPCR2,WIFI,H,IT,DE,S,1,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,IT	48047
99512-000107	CR2-2-000121	LPCR2,WIFI,H,IT,DE,S,1,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,IT	48047
99512-000108	CR2-2-000122	LPCR2,CELL,H,IT,EN-GB,S,2,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,IT	48047
99512-000109	CR2-2-000123	LPCR2,WIFI,B,FR,DE,S,1,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,FR	48047
99512-000110	CR2-2-000124	LPCR2,WIFI,H,FR,DE,A,1,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,FR	48047
99512-000111	CR2-2-000125	LPCR2,WIFI,H,FR,EN-GB,A,2,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,FR	48047
99512-000112	CR2-2-000126	LPCR2,CELL,B,FR,DE,S,1,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,FR	48047
99512-000113	CR2-2-000127	LPCR2,WIFI,H,ES-ES,DE,S,1,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,ES	48047
99512-000114	CR2-2-000128	LPCR2,WIFI,H,PL,DE,S,1,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,PL	48047
99512-000115	CR2-2-000129	LPCR2,WIFI,H,PL,DE,S,1,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,PL	48047
99512-000116	CR2-2-000130	LPCR2,CELL,H,PL,DE,S,1,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,PL	48047
99512-000117	CR2-2-000131	LPCR2,CELL,H,PL,DE,S,1,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,PL	48047
99512-000118	CR2-2-000132	LPCR2,WIFI,B,EN-GB,DE,S,1,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,AU	48047
99512-000119	CR2-2-000133	LPCR2,WIFI,H,EN-GB,DE,S,1,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,AU	48047
99512-000120	CR2-2-000134	LPCR2,WIFI,B,EN-GB,DE,A,1,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,AU	48047
99512-000121	CR2-2-000135	LPCR2,WIFI,H,EN-GB,DE,A,1,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,AU	48047
99512-000122	CR2-2-000136	LPCR2,WIFI,B,EN-GB,DE,S,1,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,NZ	48047
99512-000123	CR2-2-000137	LPCR2,WIFI,H,EN-GB,DE,S,1,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,NZ	48047
99512-000124	CR2-2-000138	LPCR2,WIFI,B,EN-GB,DE,A,1,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,NZ	48047
99512-000125	CR2-2-000139	LPCR2,WIFI,H,EN-GB,DE,A,1,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,NZ	48047
99512-000126	CR2-2-000140	LPCR2,USB,H,ES-ES,DE,S,1,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,GT	48047
99512-000127	CR2-2-000141	LPCR2,USB,H,EN-GB,DE,A,1,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,ZA	48047
99512-000128	CR2-2-000142	LPCR2,WIFI,H,EN-GB,DE,A,1,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,AE	48047
99512-000129	CR2-2-000143	LPCR2,WIFI,H,EN-GB,DE,A,1,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,QA	48047
99512-000130	CR2-2-000144	LPCR2,WIFI,H,EN-GB,DE,S,1,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,MT	48047
99512-000131	CR2-2-000145	LPCR2,USB,H,EN-GB,DE,S,1,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,MY	48047
99512-000132	CR2-2-000146	LPCR2,USB,H,EN-GB,DE,S,1,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,PH	48047
99512-000133	CR2-2-000147	LPCR2,WIFI,H,NO,EN-GB,S,2,ROS,MODET,ROS,FLX,CKB,CPR180,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,NO	48047
99512-000134	CR2-2-000148	LPCR2,CELL,B,NO,DE,S,1,ROS,MODET,ROS,FLX,CKB,CPR180,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,NO	48047
99512-000135	CR2-2-000149	LPCR2,CELL,H,NO,DE,S,1,ROS,MODET,ROS,FLX,CKB,CPR180,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,NO	48047
99512-000136	CR2-2-000150	LPCR2,CELL,B,NO,DE,A,1,ROS,MODET,ROS,FLX,CKB,CPR180,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,NO	48047
99512-000137	CR2-2-000151	LPCR2,CELL,H,NO,DE,A,1,ROS,MODET,ROS,FLX,CKB,CPR180,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,NO	48047
99512-000138	CR2-2-000152	LPCR2,CELL,B,NO,EN-GB,S,2,ROS,MODET,ROS,FLX,CKB,CPR180,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,NO	48047
99512-000139	CR2-2-000153	LPCR2,CELL,H,NO,EN-GB,S,2,ROS,MODET,ROS,FLX,CKB,CPR180,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,NO	48047
99512-000140	CR2-2-000154	LPCR2,WIFI,B,FI,DE,S,1,ROS,MODET,ROS,FLX,CKB,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,FI	48047

Ref #	Part Number	Description	GMDN
99512-000273	CR2-2-000287	LPCR2,CELL,B,IT,EN-GB,S,2,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,CH	48047
99512-000274	CR2-2-000288	LPCR2,CELL,H,IT,EN-GB,S,2,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,CH	48047
99512-000275	CR2-2-000289	LPCR2,CELL,B,FR,EN-GB,S,2,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,CH	48047
99512-000276	CR2-2-000290	LPCR2,CELL,H,FR,EN-GB,S,2,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,CH	48047
99512-000277	CR2-2-000291	LPCR2,WIFI,B,FR,DE,S,1,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,CH	48047
99512-000278	CR2-2-000292	LPCR2,CELL,B,FR,EN-GB,A,2,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,CH	48047
99512-000279	CR2-2-000293	LPCR2,WIFI,H,FR,DE,S,1,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,CH	48047
99512-000280	CR2-2-000294	LPCR2,WIFI,B,FR,DE,A,1,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,CH	48047
99512-000281	CR2-2-000295	LPCR2,CELL,H,FR,EN-GB,A,2,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,CH	48047
99512-000282	CR2-2-000296	LPCR2,WIFI,H,FR,DE,A,1,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,CH	48047
99512-000283	CR2-2-000297	LPCR2,WIFI,B,FR,DE,S,1,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,CH	48047
99512-000284	CR2-2-000298	LPCR2,WIFI,H,FR,DE,S,1,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,CH	48047
99512-000285	CR2-2-000299	LPCR2,WIFI,B,FR,DE,A,1,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,CH	48047
99512-000286	CR2-2-000300	LPCR2,WIFI,H,FR,DE,A,1,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,CH	48047
99512-000287	CR2-2-000301	LPCR2,WIFI,B,FR,EN-GB,S,2,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,CH	48047
99512-000288	CR2-2-000302	LPCR2,CELL,B,FR,EN-GB,S,2,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,CH	48047
99512-000289	CR2-2-000303	LPCR2,WIFI,H,FR,EN-GB,S,2,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,CH	48047
99512-000290	CR2-2-000304	LPCR2,WIFI,B,FR,EN-GB,S,2,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,CH	48047
99512-000291	CR2-2-000305	LPCR2,CELL,H,FR,EN-GB,S,2,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,CH	48047
99512-000292	CR2-2-000306	LPCR2,WIFI,H,FR,EN-GB,S,2,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,CH	48047
99512-000293	CR2-2-000307	LPCR2,CELL,B,FR,DE,S,1,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,CH	48047
99512-000294	CR2-2-000308	LPCR2,CELL,H,FR,DE,S,1,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,CH	48047
99512-000295	CR2-2-000309	LPCR2,WIFI,B,DE,EN-GB,S,1,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,CH	48047
99512-000296	CR2-2-000310	LPCR2,CELL,B,FR,DE,A,1,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,CH	48047
99512-000297	CR2-2-000311	LPCR2,CELL,H,FR,DE,A,1,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,CH	48047
99512-000298	CR2-2-000312	LPCR2,WIFI,H,DE,EN-GB,S,1,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,CH	48047
99512-000299	CR2-2-000313	LPCR2,WIFI,B,DE,EN-GB,A,1,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,CH	48047
99512-000300	CR2-2-000314	LPCR2,WIFI,H,DE,EN-GB,A,1,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,CH	48047
99512-000301	CR2-2-000315	LPCR2,WIFI,B,DE,EN-GB,S,1,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,CH	48047
99512-000302	CR2-2-000316	LPCR2,WIFI,H,DE,EN-GB,S,1,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,CH	48047
99512-000303	CR2-2-000317	LPCR2,WIFI,B,DE,EN-GB,A,1,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,CH	48047
99512-000304	CR2-2-000318	LPCR2,WIFI,H,DE,EN-GB,A,1,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,CH	48047
99512-000305	CR2-2-000319	LPCR2,WIFI,B,DE,EN-GB,S,2,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,CH	48047
99512-000306	CR2-2-000320	LPCR2,WIFI,H,DE,EN-GB,S,2,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,CH	48047
99512-000307	CR2-2-000321	LPCR2,WIFI,B,DE,EN-GB,S,2,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,CH	48047
99512-000308	CR2-2-000322	LPCR2,WIFI,H,DE,EN-GB,S,2,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,CH	48047
99512-000309	CR2-2-000323	LPCR2,CELL,B,DE,EN-GB,S,1,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,CH	48047
99512-000310	CR2-2-000324	LPCR2,CELL,H,IT,DE,S,1,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,CH	48047
99512-000311	CR2-2-000325	LPCR2,CELL,B,DE,EN-GB,A,1,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,CH	48047
99512-000312	CR2-2-000326	LPCR2,CELL,H,DE,EN-GB,A,1,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,CH	48047
99512-000313	CR2-2-000327	LPCR2,CELL,B,DE,EN-GB,S,1,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,CH	48047
99512-000314	CR2-2-000328	LPCR2,CELL,H,DE,EN-GB,S,1,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,CH	48047
99512-000315	CR2-2-000329	LPCR2,CELL,B,DE,EN-GB,A,1,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,CH	48047
99512-000316	CR2-2-000330	LPCR2,CELL,H,DE,EN-GB,A,1,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,CH	48047

Ref #	Part Number	Description	GMDN
99512-000405	CR2-2-000419	LPCR2,WIFI,H,DE,EN-GB,S,1,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,LI	48047
99512-000406	CR2-2-000420	LPCR2,WIFI,B,DE,EN-GB,A,1,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,LI	48047
99512-000407	CR2-2-000421	LPCR2,WIFI,H,DE,EN-GB,A,1,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,LI	48047
99512-000408	CR2-2-000422	LPCR2,WIFI,B,DE,EN-GB,S,1,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,LI	48047
99512-000409	CR2-2-000423	LPCR2,USB,H,ZH-HK,DE,A,1,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,HK	48047
99512-000410	CR2-2-000424	LPCR2,WIFI,H,DE,EN-GB,S,1,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,LI	48047
99512-000411	CR2-2-000425	LPCR2,WIFI,B,DE,EN-GB,A,1,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,LI	48047
99512-000412	CR2-2-000426	LPCR2,WIFI,H,DE,EN-GB,A,1,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,LI	48047
99512-000413	CR2-2-000427	LPCR2,WIFI,B,DE,EN-GB,S,2,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,LI	48047
99512-000414	CR2-2-000428	LPCR2,USB,B,ZH-HK,EN-GB,S,2,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,HK	48047
99512-000415	CR2-2-000429	LPCR2,WIFI,H,DE,EN-GB,S,2,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,LI	48047
99512-000416	CR2-2-000430	LPCR2,WIFI,B,DE,EN-GB,A,2,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,LI	48047
99512-000417	CR2-2-000431	LPCR2,WIFI,H,DE,EN-GB,A,2,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,LI	48047
99512-000418	CR2-2-000432	LPCR2,WIFI,B,DE,EN-GB,S,2,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,LI	48047
99512-000419	CR2-2-000433	LPCR2,WIFI,H,DE,EN-GB,S,2,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,LI	48047
99512-000420	CR2-2-000434	LPCR2,CELL,H,DE,EN-GB,S,1,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,LI	48047
99512-000421	CR2-2-000435	LPCR2,CELL,B,DE,EN-GB,S,1,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,LI	48047
99512-000422	CR2-2-000436	LPCR2,CELL,B,DE,EN-GB,A,1,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,LI	48047
99512-000423	CR2-2-000437	LPCR2,CELL,H,DE,EN-GB,A,1,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,LI	48047
99512-000424	CR2-2-000438	LPCR2,CELL,B,DE,EN-GB,S,1,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,LI	48047
99512-000425	CR2-2-000439	LPCR2,CELL,H,DE,EN-GB,S,1,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,LI	48047
99512-000426	CR2-2-000440	LPCR2,CELL,B,DE,EN-GB,A,1,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,LI	48047
99512-000427	CR2-2-000441	LPCR2,USB,H,ZH-HK,EN-GB,S,2,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,HK	48047
99512-000428	CR2-2-000442	LPCR2,CELL,H,DE,EN-GB,A,1,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,LI	48047
99512-000429	CR2-2-000443	LPCR2,CELL,B,DE,EN-GB,S,2,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,LI	48047
99512-000430	CR2-2-000444	LPCR2,CELL,H,DE,EN-GB,S,2,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,LI	48047
99512-000431	CR2-2-000445	LPCR2,CELL,B,DE,EN-GB,A,2,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,LI	48047
99512-000432	CR2-2-000446	LPCR2,USB,B,ZH-HK,EN-GB,A,2,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,HK	48047
99512-000433	CR2-2-000447	LPCR2,CELL,H,DE,EN-GB,A,2,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,LI	48047
99512-000434	CR2-2-000448	LPCR2,CELL,B,DE,EN-GB,S,2,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,LI	48047
99512-000435	CR2-2-000449	LPCR2,USB,H,ZH-HK,EN-GB,A,2,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,HK	48047
99512-000436	CR2-2-000450	LPCR2,CELL,H,DE,EN-GB,S,2,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,LI	48047
99512-000437	CR2-2-000451	LPCR2,WIFI,B,DE,EN-GB,S,1,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,AT	48047
99512-000438	CR2-2-000452	LPCR2,WIFI,H,DE,EN-GB,S,1,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,AT	48047
99512-000439	CR2-2-000453	LPCR2,WIFI,B,DE,EN-GB,S,1,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,AT	48047
99512-000440	CR2-2-000454	LPCR2,WIFI,H,DE,EN-GB,S,1,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,AT	48047
99512-000441	CR2-2-000455	LPCR2,WIFI,B,DE,EN-GB,S,2,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,AT	48047
99512-000442	CR2-2-000456	LPCR2,WIFI,H,DE,EN-GB,S,2,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,AT	48047
99512-000443	CR2-2-000457	LPCR2,WIFI,B,DE,EN-GB,S,2,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,AT	48047
99512-000444	CR2-2-000458	LPCR2,USB,B,EN-GB,DE,S,1,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,HK	48047
99512-000445	CR2-2-000459	LPCR2,WIFI,H,DE,EN-GB,S,2,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,AT	48047
99512-000446	CR2-2-000460	LPCR2,CELL,B,DE,EN-GB,S,1,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,AT	48047
99512-000447	CR2-2-000461	LPCR2,CELL,H,DE,EN-GB,S,1,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,AT	48047
99512-000448	CR2-2-000462	LPCR2,USB,H,EN-GB,DE,S,1,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,HK	48047

Ref #	Part Number	Description	GMDN
99512-000449	CR2-2-000463	LPCR2,CELL,B,DE,EN-GB,A,1,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,AT	48047
99512-000450	CR2-2-000464	LPCR2,CELL,H,DE,EN-GB,A,1,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,AT	48047
99512-000451	CR2-2-000465	LPCR2,CELL,B,DE,EN-GB,S,1,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,AT	48047
99512-000452	CR2-2-000466	LPCR2,USB,B,EN-GB,DE,A,1,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,HK	48047
99512-000453	CR2-2-000467	LPCR2,CELL,H,DE,EN-GB,S,1,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,AT	48047
99512-000454	CR2-2-000468	LPCR2,CELL,B,DE,EN-GB,A,1,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,AT	48047
99512-000455	CR2-2-000469	LPCR2,CELL,H,DE,EN-GB,A,1,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,AT	48047
99512-000456	CR2-2-000470	LPCR2,CELL,B,DE,EN-GB,S,2,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,AT	48047
99512-000457	CR2-2-000471	LPCR2,CELL,H,DE,EN-GB,S,2,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,AT	48047
99512-000458	CR2-2-000472	LPCR2,CELL,B,DE,EN-GB,A,2,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,AT	48047
99512-000459	CR2-2-000473	LPCR2,CELL,H,DE,EN-GB,A,2,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,AT	48047
99512-000460	CR2-2-000474	LPCR2,CELL,B,DE,EN-GB,S,2,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,AT	48047
99512-000461	CR2-2-000475	LPCR2,CELL,H,DE,EN-GB,S,2,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,AT	48047
99512-000462	CR2-2-000476	LPCR2,CELL,B,EN-GB,DE,S,1,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,NZ	48047
99512-000463	CR2-2-000477	LPCR2,CELL,H,EN-GB,DE,S,1,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,NZ	48047
99512-000464	CR2-2-000478	LPCR2,CELL,B,EN-GB,DE,A,1,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,NZ	48047
99512-000465	CR2-2-000479	LPCR2,CELL,H,EN-GB,DE,A,1,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,NZ	48047
99512-000466	CR2-2-000480	LPCR2,WIFI,B,EN-GB,KO,S,2,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,NZ	48047
99512-000467	CR2-2-000481	LPCR2,WIFI,H,EN-GB,KO,S,2,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,NZ	48047
99512-000468	CR2-2-000482	LPCR2,WIFI,B,EN-GB,KO,A,2,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,NZ	48047
99512-000469	CR2-2-000483	LPCR2,WIFI,H,EN-GB,KO,A,2,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,NZ	48047
99512-000470	CR2-2-000484	LPCR2,CELL,B,EN-GB,KO,S,2,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,NZ	48047
99512-000471	CR2-2-000485	LPCR2,CELL,H,EN-GB,KO,S,2,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,NZ	48047
99512-000472	CR2-2-000486	LPCR2,CELL,B,EN-GB,KO,A,2,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,NZ	48047
99512-000473	CR2-2-000487	LPCR2,CELL,H,EN-GB,KO,A,2,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,NZ	48047
99512-000474	CR2-2-000488	LPCR2,WIFI,B,IT,DE,S,1,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,IT	48047
99512-000475	CR2-2-000489	LPCR2,WIFI,B,IT,DE,A,1,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,IT	48047
99512-000476	CR2-2-000490	LPCR2,WIFI,H,IT,DE,A,1,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,IT	48047
99512-000477	CR2-2-000491	LPCR2,WIFI,B,IT,DE,S,1,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,IT	48047
99512-000478	CR2-2-000492	LPCR2,WIFI,B,IT,EN-GB,S,2,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,IT	48047
99512-000479	CR2-2-000493	LPCR2,WIFI,H,IT,EN-GB,S,2,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,IT	48047
99512-000480	CR2-2-000494	LPCR2,WIFI,B,IT,EN-GB,A,2,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,IT	48047
99512-000481	CR2-2-000495	LPCR2,WIFI,H,IT,EN-GB,A,2,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,IT	48047
99512-000482	CR2-2-000496	LPCR2,WIFI,B,IT,EN-GB,S,2,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,IT	48047
99512-000483	CR2-2-000497	LPCR2,WIFI,H,IT,EN-GB,S,2,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,IT	48047
99512-000484	CR2-2-000498	LPCR2,CELL,B,IT,DE,S,1,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,IT	48047
99512-000485	CR2-2-000499	LPCR2,CELL,H,IT,DE,S,1,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,IT	48047
99512-000486	CR2-2-000500	LPCR2,CELL,B,IT,DE,A,1,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,IT	48047
99512-000487	CR2-2-000501	LPCR2,CELL,H,IT,DE,A,1,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,IT	48047
99512-000488	CR2-2-000502	LPCR2,CELL,B,IT,DE,S,1,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,IT	48047
99512-000489	CR2-2-000503	LPCR2,CELL,H,IT,DE,S,1,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,IT	48047
99512-000490	CR2-2-000504	LPCR2,CELL,B,IT,EN-GB,S,2,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,IT	48047
99512-000491	CR2-2-000505	LPCR2,CELL,H,IT,EN-GB,S,2,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,IT	48047
99512-000492	CR2-2-000506	LPCR2,CELL,B,IT,EN-GB,A,2,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,IT	48047

Ref #	Part Number	Description	GMDN
99512-000493	CR2-2-000507	LPCR2,CELL,H,IT,EN-GB,A,2,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,IT	48047
99512-000494	CR2-2-000508	LPCR2,CELL,B,IT,EN-GB,S,2,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,IT	48047
99512-000495	CR2-2-000509	LPCR2,WIFI,H,FR,DE,S,1,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,FR	48047
99512-000496	CR2-2-000510	LPCR2,WIFI,B,FR,DE,A,1,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,FR	48047
99512-000497	CR2-2-000511	LPCR2,WIFI,B,FR,DE,S,1,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,FR	48047
99512-000498	CR2-2-000512	LPCR2,WIFI,H,FR,DE,S,1,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,FR	48047
99512-000499	CR2-2-000513	LPCR2,WIFI,B,FR,DE,A,1,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,FR	48047
99512-000500	CR2-2-000514	LPCR2,WIFI,H,FR,DE,A,1,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,FR	48047
99512-000501	CR2-2-000515	LPCR2,WIFI,B,FR,EN-GB,S,2,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,FR	48047
99512-000502	CR2-2-000516	LPCR2,WIFI,H,FR,EN-GB,S,2,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,FR	48047
99512-000503	CR2-2-000517	LPCR2,WIFI,B,FR,EN-GB,A,2,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,FR	48047
99512-000504	CR2-2-000518	LPCR2,WIFI,B,FR,EN-GB,S,2,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,FR	48047
99512-000505	CR2-2-000519	LPCR2,WIFI,H,FR,EN-GB,S,2,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,FR	48047
99512-000506	CR2-2-000520	LPCR2,CELL,H,FR,DE,S,1,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,FR	48047
99512-000507	CR2-2-000521	LPCR2,CELL,B,FR,DE,A,1,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,FR	48047
99512-000508	CR2-2-000522	LPCR2,CELL,H,FR,DE,A,1,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,FR	48047
99512-000509	CR2-2-000523	LPCR2,CELL,B,FR,DE,S,1,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,FR	48047
99512-000510	CR2-2-000524	LPCR2,CELL,H,FR,DE,S,1,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,FR	48047
99512-000511	CR2-2-000525	LPCR2,CELL,B,FR,DE,A,1,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,FR	48047
99512-000512	CR2-2-000526	LPCR2,CELL,H,FR,DE,A,1,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,FR	48047
99512-000513	CR2-2-000527	LPCR2,CELL,B,FR,EN-GB,S,2,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,FR	48047
99512-000029	CR2-2-000528	LPCR2,CELL,H,FR,EN-GB,S,ROS,MODET,ROSON,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,FR	48047
99512-000514	CR2-2-000529	LPCR2,CELL,B,FR,EN-GB,A,2,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,FR	48047
99512-000515	CR2-2-000530	LPCR2,CELL,H,FR,EN-GB,A,2,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,FR	48047
99512-000516	CR2-2-000531	LPCR2,CELL,B,FR,EN-GB,S,2,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,FR	48047
99512-000517	CR2-2-000532	LPCR2,CELL,H,FR,EN-GB,S,2,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,FR	48047
99512-000518	CR2-2-000533	LPCR2,USB,B,FR,DE,S,1,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,GF	48047
99512-000519	CR2-2-000534	LPCR2,USB,H,FR,DE,S,1,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,GF	48047
99512-000520	CR2-2-000535	LPCR2,USB,B,FR,DE,A,1,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,GF	48047
99512-000521	CR2-2-000536	LPCR2,USB,H,FR,DE,A,1,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,GF	48047
99512-000522	CR2-2-000537	LPCR2,USB,B,FR,DE,S,1,MODET,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,GF	48047
99512-000523	CR2-2-000538	LPCR2,USB,H,FR,DE,S,1,MODET,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,GF	48047
99512-000524	CR2-2-000539	LPCR2,USB,B,FR,DE,A,1,MODET,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,GF	48047
99512-000525	CR2-2-000540	LPCR2,USB,H,FR,DE,A,1,MODET,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,GF	48047
99512-000526	CR2-2-000541	LPCR2,USB,B,FR,EN-GB,S,2,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,GF	48047
99512-000527	CR2-2-000542	LPCR2,USB,H,FR,EN-GB,S,2,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,GF	48047
99512-000528	CR2-2-000543	LPCR2,USB,B,FR,EN-GB,A,2,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,GF	48047
99512-000529	CR2-2-000544	LPCR2,USB,H,FR,EN-GB,A,2,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,GF	48047
99512-000530	CR2-2-000545	LPCR2,USB,B,FR,EN-GB,S,2,MODET,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,GF	48047
99512-000531	CR2-2-000546	LPCR2,USB,H,FR,EN-GB,S,2,MODET,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,GF	48047
99512-000532	CR2-2-000547	LPCR2,USB,B,FR,DE,S,1,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,PF	48047
99512-000533	CR2-2-000548	LPCR2,USB,H,FR,DE,S,1,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,PF	48047
99512-000534	CR2-2-000549	LPCR2,USB,B,FR,DE,A,1,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,PF	48047
99512-000535	CR2-2-000550	LPCR2,USB,H,FR,DE,A,1,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,PF	48047

Ref #	Part Number	Description	GMDN
99512-000536	CR2-2-000551	LPCR2,USB,B,FR,DE,S,1,MODET,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,PF	48047
99512-000537	CR2-2-000552	LPCR2,USB,H,FR,DE,S,1,MODET,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,PF	48047
99512-000538	CR2-2-000553	LPCR2,USB,B,FR,DE,A,1,MODET,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,PF	48047
99512-000539	CR2-2-000554	LPCR2,USB,H,FR,DE,A,1,MODET,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,PF	48047
99512-000540	CR2-2-000555	LPCR2,USB,B,FR,EN-GB,S,2,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,PF	48047
99512-000541	CR2-2-000556	LPCR2,USB,H,FR,EN-GB,S,2,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,PF	48047
99512-000542	CR2-2-000557	LPCR2,USB,B,FR,EN-GB,A,2,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,PF	48047
99512-000543	CR2-2-000558	LPCR2,USB,H,FR,EN-GB,A,2,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,PF	48047
99512-000544	CR2-2-000559	LPCR2,USB,B,FR,EN-GB,S,2,MODET,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,PF	48047
99512-000545	CR2-2-000560	LPCR2,USB,H,FR,EN-GB,S,2,MODET,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,PF	48047
99512-000546	CR2-2-000561	LPCR2,WIFI,B,ES-ES,DE,S,1,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,ES	48047
99512-000547	CR2-2-000562	LPCR2,WIFI,B,ES-ES,EN-GB,S,2,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,ES	48047
99512-000548	CR2-2-000563	LPCR2,WIFI,H,ES-ES,EN-GB,S,2,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,ES	48047
99512-000549	CR2-2-000564	LPCR2,WIFI,B,ES-ES,EN-GB,A,2,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,ES	48047
99512-000550	CR2-2-000565	LPCR2,WIFI,H,ES-ES,EN-GB,A,2,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,ES	48047
99512-000551	CR2-2-000566	LPCR2,CELL,B,ES-ES,DE,S,1,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,ES	48047
99512-000552	CR2-2-000567	LPCR2,CELL,H,ES-ES,DE,S,1,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,ES	48047
99512-000553	CR2-2-000568	LPCR2,CELL,B,ES-ES,EN-GB,S,2,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,ES	48047
99512-000554	CR2-2-000569	LPCR2,CELL,H,ES-ES,EN-GB,S,2,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,ES	48047
99512-000555	CR2-2-000570	LPCR2,CELL,B,ES-ES,EN-GB,A,2,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,ES	48047
99512-000556	CR2-2-000571	LPCR2,CELL,H,ES-ES,EN-GB,A,2,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,ES	48047
99512-000557	CR2-2-000572	LPCR2,WIFI,B,PL,DE,S,1,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,PL	48047
99512-000558	CR2-2-000573	LPCR2,WIFI,B,PL,DE,A,1,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,PL	48047
99512-000559	CR2-2-000574	LPCR2,WIFI,H,PL,DE,A,1,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,PL	48047
99512-000560	CR2-2-000575	LPCR2,WIFI,B,PL,DE,S,1,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,PL	48047
99512-000561	CR2-2-000576	LPCR2,WIFI,B,PL,DE,A,1,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,PL	48047
99512-000562	CR2-2-000577	LPCR2,WIFI,H,PL,DE,A,1,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,PL	48047
99512-000563	CR2-2-000578	LPCR2,WIFI,B,PL,EN-GB,S,2,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,PL	48047
99512-000564	CR2-2-000579	LPCR2,WIFI,H,PL,EN-GB,S,2,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,PL	48047
99512-000565	CR2-2-000580	LPCR2,WIFI,B,PL,EN-GB,S,2,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,PL	48047
99512-000566	CR2-2-000581	LPCR2,WIFI,H,PL,EN-GB,S,2,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,PL	48047
99512-000567	CR2-2-000582	LPCR2,WIFI,B,PL,DE,S,2,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,PL	48047
99512-000568	CR2-2-000583	LPCR2,WIFI,H,PL,DE,S,2,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,PL	48047
99512-000569	CR2-2-000584	LPCR2,WIFI,B,PL,DE,S,2,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,PL	48047
99512-000570	CR2-2-000585	LPCR2,WIFI,H,PL,DE,S,2,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,PL	48047
99512-000571	CR2-2-000586	LPCR2,USB,H,EN-GB,DE,A,1,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,HK	48047
99512-000572	CR2-2-000587	LPCR2,CELL,B,PL,DE,S,1,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,PL	48047
99512-000573	CR2-2-000588	LPCR2,CELL,B,PL,DE,A,1,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,PL	48047
99512-000574	CR2-2-000589	LPCR2,CELL,H,PL,DE,A,1,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,PL	48047
99512-000575	CR2-2-000590	LPCR2,CELL,B,PL,DE,S,1,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,PL	48047
99512-000576	CR2-2-000591	LPCR2,CELL,B,PL,DE,A,1,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,PL	48047
99512-000577	CR2-2-000592	LPCR2,CELL,H,PL,DE,A,1,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,PL	48047
99512-000578	CR2-2-000593	LPCR2,CELL,B,PL,EN-GB,S,2,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,PL	48047
99512-000579	CR2-2-000594	LPCR2,CELL,H,PL,EN-GB,S,2,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,PL	48047

Ref #	Part Number	Description	GMDN
99512-000580	CR2-2-000595	LPCR2,CELL,B,PL,EN-GB,S,2,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,PL	48047
99512-000581	CR2-2-000596	LPCR2,CELL,H,PL,EN-GB,S,2,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,PL	48047
99512-000582	CR2-2-000597	LPCR2,CELL,B,PL,DE,S,2,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,PL	48047
99512-000583	CR2-2-000598	LPCR2,CELL,H,PL,DE,S,2,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,PL	48047
99512-000584	CR2-2-000599	LPCR2,CELL,B,PL,DE,S,2,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,PL	48047
99512-000585	CR2-2-000600	LPCR2,USB,B,EN-GB,ZH-HK,S,2,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,HK	48047
99512-000586	CR2-2-000601	LPCR2,CELL,H,PL,DE,S,2,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,PL	48047
99512-000587	CR2-2-000602	LPCR2,USB,H,EN-GB,ZH-HK,S,2,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,HK	48047
99512-000588	CR2-2-000603	LPCR2,CELL,B,EN-GB,DE,S,1,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,AU	48047
99512-000589	CR2-2-000604	LPCR2,CELL,H,EN-GB,DE,S,1,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,AU	48047
99512-000590	CR2-2-000605	LPCR2,CELL,B,EN-GB,DE,A,1,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,AU	48047
99512-000591	CR2-2-000606	LPCR2,USB,B,EN-GB,ZH-HK,A,2,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,HK	48047
99512-000592	CR2-2-000607	LPCR2,CELL,H,EN-GB,DE,A,1,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,AU	48047
99512-000593	CR2-2-000608	LPCR2,WIFI,B,EN-GB,IT,S,2,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,AU	48047
99512-000594	CR2-2-000609	LPCR2,WIFI,H,EN-GB,IT,S,2,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,AU	48047
99512-000595	CR2-2-000610	LPCR2,WIFI,B,EN-GB,IT,A,2,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,AU	48047
99512-000596	CR2-2-000611	LPCR2,USB,H,EN-GB,ZH-HK,A,2,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,HK	48047
99512-000597	CR2-2-000612	LPCR2,WIFI,H,EN-GB,IT,A,2,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,AU	48047
99512-000598	CR2-2-000613	LPCR2,CELL,B,EN-GB,IT,S,2,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,AU	48047
99512-000599	CR2-2-000614	LPCR2,CELL,H,EN-GB,IT,S,2,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,AU	48047
99512-000600	CR2-2-000615	LPCR2,CELL,B,EN-GB,IT,A,2,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,AU	48047
99512-000601	CR2-2-000616	LPCR2,CELL,H,EN-GB,IT,A,2,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,AU	48047
99512-000602	CR2-2-000617	LPCR2,USB,B,ES-ES,DE,S,1,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,GT	48047
99512-000603	CR2-2-000618	LPCR2,USB,B,ES-ES,DE,A,1,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,GT	48047
99512-000604	CR2-2-000619	LPCR2,USB,H,ES-ES,DE,A,1,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,GT	48047
99512-000605	CR2-2-000620	LPCR2,USB,B,ES-ES,EN-GB,S,2,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,GT	48047
99512-000606	CR2-2-000621	LPCR2,USB,H,ES-ES,EN-GB,S,2,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,GT	48047
99512-000607	CR2-2-000622	LPCR2,USB,B,ES-ES,EN-GB,A,2,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,GT	48047
99512-000608	CR2-2-000623	LPCR2,USB,H,ES-ES,EN-GB,A,2,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,GT	48047
99512-000609	CR2-2-000624	LPCR2,USB,B,EN-GB,DE,S,1,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,ZA	48047
99512-000610	CR2-2-000625	LPCR2,USB,H,EN-GB,DE,S,1,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,ZA	48047
99512-000611	CR2-2-000626	LPCR2,USB,B,EN-GB,DE,A,1,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,ZA	48047
99512-000612	CR2-2-000627	LPCR2,WIFI,B,EN-GB,DE,S,1,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,AE	48047
99512-000613	CR2-2-000628	LPCR2,WIFI,H,EN-GB,DE,S,1,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,AE	48047
99512-000614	CR2-2-000629	LPCR2,WIFI,B,EN-GB,DE,A,1,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,AE	48047
99512-000615	CR2-2-000630	LPCR2,CELL,B,EN-GB,DE,S,1,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,AE	48047
99512-000616	CR2-2-000631	LPCR2,CELL,H,EN-GB,DE,S,1,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,AE	48047
99512-000617	CR2-2-000632	LPCR2,CELL,B,EN-GB,DE,A,1,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,AE	48047
99512-000618	CR2-2-000633	LPCR2,CELL,H,EN-GB,DE,A,1,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,AE	48047
99512-000619	CR2-2-000634	LPCR2,WIFI,B,EN-GB,DE,S,1,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,QA	48047
99512-000620	CR2-2-000635	LPCR2,WIFI,H,EN-GB,DE,S,1,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,QA	48047
99512-000621	CR2-2-000636	LPCR2,WIFI,B,EN-GB,DE,A,1,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,QA	48047
99512-000622	CR2-2-000637	LPCR2,CELL,B,EN-GB,DE,S,1,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,QA	48047
99512-000623	CR2-2-000638	LPCR2,CELL,H,EN-GB,DE,S,1,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,QA	48047

Ref #	Part Number	Description	GMDN
99512-000624	CR2-2-000639	LPCR2,CELL,B,EN-GB,DE,A,1,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,QA	48047
99512-000625	CR2-2-000640	LPCR2,CELL,H,EN-GB,DE,A,1,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,QA	48047
99512-000626	CR2-2-000641	LPCR2,WIFI,B,EN-GB,DE,S,1,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,MT	48047
99512-000627	CR2-2-000642	LPCR2,WIFI,B,EN-GB,DE,A,1,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,MT	48047
99512-000628	CR2-2-000643	LPCR2,WIFI,H,EN-GB,DE,A,1,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,MT	48047
99512-000629	CR2-2-000644	LPCR2,CELL,B,EN-GB,DE,S,1,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,MT	48047
99512-000630	CR2-2-000645	LPCR2,CELL,H,EN-GB,DE,S,1,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,MT	48047
99512-000631	CR2-2-000646	LPCR2,CELL,B,EN-GB,DE,A,1,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,MT	48047
99512-000632	CR2-2-000647	LPCR2,CELL,H,EN-GB,DE,A,1,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,MT	48047
99512-000633	CR2-2-000648	LPCR2,USB,B,EN-GB,DE,S,1,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,MY	48047
99512-000634	CR2-2-000649	LPCR2,USB,B,EN-GB,DE,A,1,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,MY	48047
99512-000635	CR2-2-000650	LPCR2,USB,H,EN-GB,DE,A,1,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,MY	48047
99512-000636	CR2-2-000651	LPCR2,USB,B,EN-GB,DE,S,1,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,PH	48047
99512-000637	CR2-2-000652	LPCR2,USB,B,EN-GB,DE,A,1,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,PH	48047
99512-000638	CR2-2-000653	LPCR2,USB,H,EN-GB,DE,A,1,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,PH	48047
99512-000639	CR2-2-000654	LPCR2,WIFI,B,EN-US,DE,S,1,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,CA	48047
99512-000640	CR2-2-000655	LPCR2,USB,B,ES-ES,DE,S,1,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,AR	48047
99512-000641	CR2-2-000656	LPCR2,USB,H,ES-ES,DE,S,1,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,AR	48047
99512-000642	CR2-2-000657	LPCR2,USB,B,ES-ES,DE,A,1,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,AR	48047
99512-000643	CR2-2-000658	LPCR2,USB,H,ES-ES,DE,A,1,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,AR	48047
99512-000644	CR2-2-000659	LPCR2,USB,B,ES-ES,EN-GB,S,2,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,AR	48047
99512-000645	CR2-2-000660	LPCR2,USB,H,ES-ES,EN-GB,S,2,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,AR	48047
99512-000646	CR2-2-000661	LPCR2,USB,B,ES-ES,EN-GB,A,2,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,AR	48047
99512-000647	CR2-2-000662	LPCR2,USB,H,ES-ES,EN-GB,A,2,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,AR	48047
99512-000648	CR2-2-000663	LPCR2,USB,B,ES-ES,DE,S,1,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,MX	48047
99512-000649	CR2-2-000664	LPCR2,USB,H,ES-ES,DE,S,1,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,MX	48047
99512-000650	CR2-2-000665	LPCR2,USB,B,ES-ES,DE,A,1,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,MX	48047
99512-000651	CR2-2-000666	LPCR2,USB,H,ES-ES,DE,A,1,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,MX	48047
99512-000652	CR2-2-000667	LPCR2,USB,B,ES-ES,EN-GB,S,2,MODET,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,MX	48047
99512-000653	CR2-2-000668	LPCR2,USB,H,ES-ES,EN-GB,S,2,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,MX	48047
99512-000654	CR2-2-000669	LPCR2,USB,B,ES-ES,EN-GB,A,2,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,MX	48047
99512-000655	CR2-2-000670	LPCR2,USB,H,ES-ES,EN-GB,A,2,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,MX	48047
99512-000656	CR2-2-000671	LPCR2,USB,B,ES-ES,DE,S,1,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,PE	48047
99512-000657	CR2-2-000672	LPCR2,USB,H,ES-ES,DE,S,1,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,PE	48047
99512-000658	CR2-2-000673	LPCR2,USB,B,ES-ES,DE,A,1,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,PE	48047
99512-000659	CR2-2-000674	LPCR2,USB,H,ES-ES,DE,A,1,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,PE	48047
99512-000660	CR2-2-000675	LPCR2,USB,B,ES-ES,EN-GB,S,2,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,PE	48047
99512-000661	CR2-2-000676	LPCR2,USB,H,ES-ES,EN-GB,S,2,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,PE	48047
99512-000662	CR2-2-000677	LPCR2,USB,B,ES-ES,DE,S,1,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,EC	48047
99512-000663	CR2-2-000678	LPCR2,USB,B,ES-ES,EN-GB,A,2,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,PE	48047
99512-000664	CR2-2-000679	LPCR2,USB,H,ES-ES,DE,S,1,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,EC	48047
99512-000665	CR2-2-000680	LPCR2,USB,H,ES-ES,EN-GB,A,2,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,PE	48047
99512-000666	CR2-2-000681	LPCR2,USB,B,ES-ES,DE,A,1,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,EC	48047
99512-000667	CR2-2-000682	LPCR2,USB,B,ES-ES,DE,S,1,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,CO	48047

Ref #	Part Number	Description	GMDN
99512-000668	CR2-2-000683	LPCR2,USB,H,ES-ES,DE,A,1,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,EC	48047
99512-000669	CR2-2-000684	LPCR2,USB,H,ES-ES,DE,S,1,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,CO	48047
99512-000670	CR2-2-000685	LPCR2,USB,B,ES-ES,DE,A,1,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,CO	48047
99512-000671	CR2-2-000686	LPCR2,USB,H,ES-ES,DE,A,1,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,CO	48047
99512-000672	CR2-2-000687	LPCR2,USB,B,ES-ES,EN-GB,S,2,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,EC	48047
99512-000673	CR2-2-000688	LPCR2,USB,B,ES-ES,EN-GB,S,2,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,CO	48047
99512-000674	CR2-2-000689	LPCR2,USB,H,ES-ES,EN-GB,S,2,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,CO	48047
99512-000675	CR2-2-000690	LPCR2,USB,H,ES-ES,EN-GB,S,2,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,EC	48047
99512-000676	CR2-2-000691	LPCR2,USB,B,ES-ES,EN-GB,A,2,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,EC	48047
99512-000677	CR2-2-000692	LPCR2,USB,H,ES-ES,EN-GB,A,2,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,EC	48047
99512-000678	CR2-2-000693	LPCR2,USB,B,ES-ES,EN-GB,A,2,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,CO	48047
99512-000679	CR2-2-000694	LPCR2,USB,H,ES-ES,EN-GB,A,2,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,CO	48047
99512-000680	CR2-2-000695	LPCR2,USB,B,ES-ES,DE,S,1,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,CR	48047
99512-000681	CR2-2-000696	LPCR2,USB,B,ES-ES,DE,S,1,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,VE	48047
99512-000682	CR2-2-000697	LPCR2,USB,H,ES-ES,DE,S,1,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,CR	48047
99512-000683	CR2-2-000698	LPCR2,USB,B,ES-ES,DE,A,1,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,CR	48047
99512-000684	CR2-2-000699	LPCR2,USB,H,ES-ES,DE,A,1,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,CR	48047
99512-000685	CR2-2-000700	LPCR2,USB,H,ES-ES,DE,S,1,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,VE	48047
99512-000686	CR2-2-000701	LPCR2,USB,B,ES-ES,EN-GB,S,2,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,CR	48047
99512-000687	CR2-2-000702	LPCR2,USB,H,ES-ES,EN-GB,S,2,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,CR	48047
99512-000688	CR2-2-000703	LPCR2,USB,B,ES-ES,DE,A,1,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,VE	48047
99512-000689	CR2-2-000704	LPCR2,USB,B,ES-ES,EN-GB,A,2,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,CR	48047
99512-000690	CR2-2-000705	LPCR2,USB,H,ES-ES,DE,A,1,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,VE	48047
99512-000691	CR2-2-000706	LPCR2,USB,H,ES-ES,EN-GB,A,2,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,CR	48047
99512-000692	CR2-2-000707	LPCR2,USB,B,ES-ES,EN-GB,S,2,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,VE	48047
99512-000693	CR2-2-000708	LPCR2,USB,H,ES-ES,EN-GB,S,2,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,VE	48047
99512-000694	CR2-2-000709	LPCR2,USB,B,ES-ES,EN-GB,A,2,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,VE	48047
99512-000695	CR2-2-000710	LPCR2,USB,H,ES-ES,EN-GB,A,2,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,VE	48047
99512-000696	CR2-2-000711	LPCR2,USB,B,KO,DE,S,1,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,KR	48047
99512-000697	CR2-2-000712	LPCR2,USB,H,KO,DE,S,1,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,KR	48047
99512-000698	CR2-2-000713	LPCR2,USB,B,KO,DE,A,1,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,KR	48047
99512-000699	CR2-2-000714	LPCR2,USB,H,KO,DE,A,1,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,KR	48047
99512-000700	CR2-2-000715	LPCR2,USB,B,KO,EN-GB,S,2,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,KR	48047
99512-000701	CR2-2-000716	LPCR2,WIFI,B,FR,DE,S,1,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,CA-Q	48047
99512-000702	CR2-2-000717	LPCR2,USB,H,KO,EN-GB,S,2,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,KR	48047
99512-000703	CR2-2-000718	LPCR2,WIFI,H,FR,DE,A,1,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,CA-Q	48047
99512-000704	CR2-2-000719	LPCR2,WIFI,B,FR,DE,S,1,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,CA-Q	48047
99512-000705	CR2-2-000720	LPCR2,WIFI,H,FR,DE,S,1,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,CA-Q	48047
99512-000706	CR2-2-000721	LPCR2,WIFI,B,FR,EN-US,S,2,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,CA-Q	48047
99512-000707	CR2-2-000722	LPCR2,WIFI,H,FR,EN-US,S,2,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,CA-Q	48047
99512-000708	CR2-2-000723	LPCR2,WIFI,B,FR,EN-US,A,2,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,CA-Q	48047
99512-000709	CR2-2-000724	LPCR2,WIFI,H,FR,EN-US,A,2,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,CA-Q	48047
99512-000710	CR2-2-000725	LPCR2,WIFI,B,FR,EN-US,S,2,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,CA-Q	48047
99512-000711	CR2-2-000726	LPCR2,WIFI,H,FR,EN-US,S,2,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,CA-Q	48047

Ref #	Part Number	Description	GMDN
99512-000712	CR2-2-000727	LPCR2,CELL,B,FR,DE,S,1,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,CA-Q	48047
99512-000713	CR2-2-000728	LPCR2,CELL,H,FR,DE,S,1,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,CA-Q	48047
99512-000714	CR2-2-000729	LPCR2,CELL,B,FR,DE,A,1,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,CA-Q	48047
99512-000715	CR2-2-000730	LPCR2,CELL,H,FR,DE,A,1,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,CA-Q	48047
99512-000716	CR2-2-000731	LPCR2,CELL,B,FR,DE,S,1,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,CA-Q	48047
99512-000717	CR2-2-000732	LPCR2,CELL,H,FR,DE,S,1,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,CA-Q	48047
99512-000718	CR2-2-000733	LPCR2,CELL,B,FR,EN-US,S,2,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,CA-Q	48047
99512-000719	CR2-2-000734	LPCR2,CELL,H,FR,EN-US,S,2,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,CA-Q	48047
99512-000720	CR2-2-000735	LPCR2,CELL,B,FR,EN-US,A,2,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,CA-Q	48047
99512-000721	CR2-2-000736	LPCR2,CELL,H,FR,EN-US,A,2,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,CA-Q	48047
99512-000722	CR2-2-000737	LPCR2,WIFI,H,EN-US,DE,S,1,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,CA	48047
99512-000723	CR2-2-000738	LPCR2,CELL,B,FR,EN-US,S,2,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,CA-Q	48047
99512-000724	CR2-2-000739	LPCR2,CELL,H,FR,EN-US,S,2,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,CA-Q	48047
99512-000725	CR2-2-000740	LPCR2,WIFI,B,EN-US,DE,A,1,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,CA	48047
99512-000726	CR2-2-000741	LPCR2,WIFI,H,EN-US,DE,A,1,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,CA	48047
99512-000727	CR2-2-000742	LPCR2,WIFI,B,EN-US,DE,S,1,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,CA	48047
99512-000728	CR2-2-000743	LPCR2,WIFI,H,EN-US,DE,S,1,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,CA	48047
99512-000729	CR2-2-000744	LPCR2,WIFI,B,EN-US,FR,S,2,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,CA	48047
99512-000730	CR2-2-000745	LPCR2,WIFI,H,EN-US,FR,S,2,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,CA	48047
99512-000731	CR2-2-000746	LPCR2,WIFI,B,EN-US,FR,A,2,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,CA	48047
99512-000732	CR2-2-000747	LPCR2,WIFI,H,EN-US,FR,A,2,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,CA	48047
99512-000733	CR2-2-000748	LPCR2,WIFI,B,EN-US,FR,S,2,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,CA	48047
99512-000734	CR2-2-000749	LPCR2,WIFI,H,EN-US,FR,S,2,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,CA	48047
99512-000735	CR2-2-000750	LPCR2,CELL,B,EN-US,DE,S,1,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,CA	48047
99512-000736	CR2-2-000751	LPCR2,CELL,H,EN-US,DE,S,1,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,CA	48047
99512-000737	CR2-2-000752	LPCR2,USB,B,KO,EN-GB,A,2,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,KR	48047
99512-000738	CR2-2-000753	LPCR2,CELL,B,EN-US,DE,A,1,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,CA	48047
99512-000739	CR2-2-000754	LPCR2,USB,H,EN-GB,DE,A,1,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,KZ	48047
99512-000740	CR2-2-000755	LPCR2,CELL,H,EN-US,DE,A,1,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,CA	48047
99512-000741	CR2-2-000756	LPCR2,USB,H,KO,EN-GB,A,2,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,KR	48047
99512-000742	CR2-2-000757	LPCR2,CELL,B,EN-US,DE,S,1,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,CA	48047
99512-000743	CR2-2-000758	LPCR2,CELL,H,EN-US,DE,S,1,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,CA	48047
99512-000744	CR2-2-000759	LPCR2,CELL,B,EN-US,FR,S,2,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,CA	48047
99512-000027	CR2-2-000760	LPCR2,CELL,H,EN-US,FR,S,ROS,MODET,ROSON,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,CA	48047
99512-000745	CR2-2-000761	LPCR2,USB,B,EN-GB,DE,A,1,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,KZ	48047
99512-000746	CR2-2-000762	LPCR2,CELL,B,EN-US,FR,A,2,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,CA	48047
99512-000747	CR2-2-000763	LPCR2,CELL,H,EN-US,FR,A,2,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,CA	48047
99512-000748	CR2-2-000764	LPCR2,CELL,B,EN-US,FR,S,2,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,CA	48047
99512-000749	CR2-2-000765	LPCR2,USB,H,EN-GB,DE,S,1,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,KZ	48047
99512-000750	CR2-2-000766	LPCR2,CELL,H,EN-US,FR,S,2,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,CA	48047
99512-000701	CR2-2-000767	LPCR2,WIFI,B,FR,DE,S,1,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,CA-Q	48047
99512-000751	CR2-2-000768	LPCR2,WIFI,H,FR,DE,S,1,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,CA-Q	48047
99512-000752	CR2-2-000769	LPCR2,WIFI,B,EN-GB,DE,S,1,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,SG	48047
99512-000753	CR2-2-000770	LPCR2,WIFI,H,EN-GB,DE,S,1,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,SG	48047

Ref #	Part Number	Description	GMDN
99512-000754	CR2-2-000771	LPCR2,USB,B,EN-GB,DE,S,1,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,KZ	48047
99512-000755	CR2-2-000772	LPCR2,WIFI,B,EN-GB,DE,A,1,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,SG	48047
99512-000756	CR2-2-000773	LPCR2,USB,H,EN-GB,DE,A,1,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,LK	48047
99512-000757	CR2-2-000774	LPCR2,WIFI,H,EN-GB,DE,A,1,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,SG	48047
99512-000758	CR2-2-000775	LPCR2,USB,B,EN-GB,DE,A,1,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,LK	48047
99512-000759	CR2-2-000776	LPCR2,CELL,B,EN-GB,DE,S,1,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,SG	48047
99512-000760	CR2-2-000777	LPCR2,USB,H,EN-GB,DE,S,1,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,LK	48047
99512-000761	CR2-2-000778	LPCR2,CELL,H,EN-GB,DE,S,1,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,SG	48047
99512-000762	CR2-2-000779	LPCR2,USB,B,EN-GB,DE,S,1,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,LK	48047
99512-000763	CR2-2-000780	LPCR2,CELL,B,EN-GB,DE,A,1,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,SG	48047
99512-000764	CR2-2-000781	LPCR2,USB,H,EN-GB,DE,A,1,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,BH	48047
99512-000765	CR2-2-000782	LPCR2,CELL,H,EN-GB,DE,A,1,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,SG	48047
99512-000766	CR2-2-000783	LPCR2,WIFI,B,EN-GB,ZH-HK,S,2,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,SG	48047
99512-000767	CR2-2-000784	LPCR2,WIFI,H,EN-GB,ZH-HK,S,2,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,SG	48047
99512-000768	CR2-2-000785	LPCR2,USB,B,EN-GB,DE,A,1,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,BH	48047
99512-000769	CR2-2-000786	LPCR2,WIFI,B,EN-GB,ZH-HK,A,2,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,SG	48047
99512-000770	CR2-2-000787	LPCR2,WIFI,H,EN-GB,ZH-HK,A,2,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,SG	48047
99512-000771	CR2-2-000788	LPCR2,USB,H,EN-GB,DE,S,1,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,BH	48047
99512-000772	CR2-2-000789	LPCR2,USB,B,EN-GB,DE,S,1,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,BH	48047
99512-000773	CR2-2-000790	LPCR2,CELL,B,EN-GB,ZH-HK,S,2,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,SG	48047
99512-000774	CR2-2-000791	LPCR2,USB,H,ES-ES,EN-GB,A,2,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,UY	48047
99512-000775	CR2-2-000792	LPCR2,CELL,H,EN-GB,ZH-HK,S,2,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,SG	48047
99512-000776	CR2-2-000793	LPCR2,USB,B,ES-ES,EN-GB,A,2,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,UY	48047
99512-000777	CR2-2-000794	LPCR2,CELL,B,EN-GB,ZH-HK,A,2,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,SG	48047
99512-000778	CR2-2-000795	LPCR2,USB,H,ES-ES,EN-GB,S,2,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,UY	48047
99512-000779	CR2-2-000796	LPCR2,CELL,H,EN-GB,ZH-HK,A,2,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,SG	48047
99512-000780	CR2-2-000797	LPCR2,USB,B,ES-ES,EN-GB,S,2,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,UY	48047
99512-000781	CR2-2-000798	LPCR2,USB,H,ES-ES,DE,A,1,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,UY	48047
99512-000782	CR2-2-000799	LPCR2,USB,B,ES-ES,DE,A,1,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,UY	48047
99512-000783	CR2-2-000800	LPCR2,USB,H,ES-ES,DE,S,1,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,UY	48047
99512-000784	CR2-2-000801	LPCR2,USB,B,ES-ES,DE,S,1,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,UY	48047
99512-000785	CR2-2-000802	LPCR2,USB,B,EN-GB,DE,S,1,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,KW	48047
99512-000786	CR2-2-000803	LPCR2,USB,H,EN-GB,DE,S,1,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,KW	48047
99512-000787	CR2-2-000804	LPCR2,USB,B,EN-GB,DE,A,1,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,KW	48047
99512-000788	CR2-2-000805	LPCR2,USB,H,EN-GB,DE,A,1,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,KW	48047
99512-000789	CR2-2-000806	LPCR2,WIFI,B,EN-GB,DE,S,1,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,IL	48047
99512-000790	CR2-2-000807	LPCR2,WIFI,H,EN-GB,DE,S,1,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,IL	48047
99512-000791	CR2-2-000808	LPCR2,WIFI,B,EN-GB,DE,A,1,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,IL	48047
99512-000792	CR2-2-000809	LPCR2,WIFI,H,EN-GB,DE,A,1,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,IL	48047
99512-000793	CR2-2-000810	LPCR2,CELL,B,EN-GB,DE,S,1,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,IL	48047
99512-000794	CR2-2-000811	LPCR2,CELL,H,EN-GB,DE,S,1,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,IL	48047
99512-000795	CR2-2-000812	LPCR2,CELL,B,EN-GB,DE,A,1,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,IL	48047
99512-000796	CR2-2-000813	LPCR2,CELL,H,EN-GB,DE,A,1,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,IL	48047
99512-000797	CR2-2-000814	LPCR2,USB,B,EN-GB,DE,S,1,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,MK	48047

Ref #	Part Number	Description	GMDN
99512-000798	CR2-2-000815	LPCR2,USB,H,EN-GB,DE,S,1,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,MK	48047
99512-000799	CR2-2-000816	LPCR2,USB,B,EN-GB,DE,A,1,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,MK	48047
99512-000800	CR2-2-000817	LPCR2,USB,H,EN-GB,DE,A,1,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,MK	48047
99512-000801	CR2-2-000818	LPCR2,USB,B,EN-GB,DE,S,1,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,ID	48047
99512-000802	CR2-2-000819	LPCR2,USB,H,EN-GB,DE,S,1,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,ID	48047
99512-000803	CR2-2-000820	LPCR2,USB,B,EN-GB,DE,A,1,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,ID	48047
99512-000804	CR2-2-000821	LPCR2,USB,B,EN-GB,DE,S,1,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,EG	48047
99512-000805	CR2-2-000822	LPCR2,USB,H,EN-GB,DE,S,1,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,EG	48047
99512-000806	CR2-2-000823	LPCR2,USB,B,EN-GB,DE,A,1,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,EG	48047
99512-000807	CR2-2-000824	LPCR2,USB,H,EN-GB,DE,A,1,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,EG	48047
99512-000808	CR2-2-000825	LPCR2,USB,H,EN-GB,DE,A,1,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,ID	48047
99512-000809	CR2-2-000826	LPCR2,WIFI,B,FR,DE,A,1,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,CA-Q	48047
99512-000810	CR2-2-000827	LPCR2,WIFI,B,PL,EN-GB,A,2,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,PL	48047
99512-000811	CR2-2-000828	LPCR2,WIFI,H,PL,EN-GB,A,2,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,PL	48047
99512-000812	CR2-2-000829	LPCR2,WIFI,B,PL,EN-GB,A,2,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,PL	48047
99512-000813	CR2-2-000830	LPCR2,WIFI,H,PL,EN-GB,A,2,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,PL	48047
99512-000814	CR2-2-000831	LPCR2,WIFI,B,PL,DE,A,2,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,PL	48047
99512-000815	CR2-2-000832	LPCR2,WIFI,H,PL,DE,A,2,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,PL	48047
99512-000816	CR2-2-000833	LPCR2,WIFI,B,PL,DE,A,2,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,PL	48047
99512-000817	CR2-2-000834	LPCR2,WIFI,H,PL,DE,A,2,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,PL	48047
99512-000818	CR2-2-000835	LPCR2,CELL,B,PL,EN-GB,A,2,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,PL	48047
99512-000819	CR2-2-000836	LPCR2,CELL,H,PL,EN-GB,A,2,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,PL	48047
99512-000820	CR2-2-000837	LPCR2,CELL,B,PL,EN-GB,A,2,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,PL	48047
99512-000821	CR2-2-000838	LPCR2,CELL,H,PL,EN-GB,A,2,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,PL	48047
99512-000822	CR2-2-000839	LPCR2,CELL,B,PL,DE,A,2,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,PL	48047
99512-000823	CR2-2-000840	LPCR2,CELL,H,PL,DE,A,2,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,PL	48047
99512-000824	CR2-2-000841	LPCR2,CELL,B,PL,DE,A,2,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,PL	48047
99512-000825	CR2-2-000842	LPCR2,CELL,H,PL,DE,A,2,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,PL	48047
99512-000826	CR2-2-000843	LPCR2,CELL,H,EN-US,FR,A,2,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,US	48047
99512-000827	CR2-2-000844	LPCR2,CELL,H,EN-GB,FR,S,2,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,GB	48047
99512-000826	CR2-2-000845	LPCR2,CELL,H,EN-US,FR,A,2,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,US	48047
99512-000018	CR2-2-000846	LPCR2,CELL,H,EN-US,FR,S,ROS,MODET,ROSON,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,US	48047
99512-000829	CR2-2-000847	LPCR2,USB,B,FR,DE,S,1,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,DZ	48047
99512-000830	CR2-2-000848	LPCR2,USB,H,FR,DE,S,1,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,DZ	48047
99512-000831	CR2-2-000849	LPCR2,USB,B,FR,DE,A,1,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,DZ	48047
99512-000832	CR2-2-000850	LPCR2,USB,H,FR,DE,A,1,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,DZ	48047
99512-000833	CR2-2-000851	LPCR2,USB,H,EN-GB,DE,S,1,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,BS	48047
99512-000834	CR2-2-000852	LPCR2,USB,B,EN-GB,DE,S,1,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,BS	48047
99512-000835	CR2-2-000853	LPCR2,USB,B,EN-GB,DE,A,1,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,BS	48047
99512-000836	CR2-2-000854	LPCR2,USB,H,EN-GB,DE,A,1,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,BS	48047
99512-000837	CR2-2-000855	LPCR2,USB,B,ES-ES,DE,S,1,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,BO	48047
99512-000838	CR2-2-000856	LPCR2,USB,H,ES-ES,DE,S,1,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,BO	48047
99512-000839	CR2-2-000857	LPCR2,USB,B,ES-ES,DE,A,1,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,BO	48047
99512-000840	CR2-2-000858	LPCR2,USB,H,ES-ES,DE,A,1,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,BO	48047

Ref #	Part Number	Description	GMDN
99512-000929	CR2-2-000947	LPCR2,CELL,H,EN-GB,NL,A,2,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,GB	48047
99512-000991	CR2-2-000948	LPCR2,WIFI,H,EN-GB,DE,S,1,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,GB	48047
99512-000930	CR2-2-000949	LPCR2,CELL,B,NL,EN-GB,S,2,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,360J,360J,50J,75J,90J,ALRT,NL	48047
99512-000931	CR2-2-000950	LPCR2,WIFI,B,EN-GB,DE,S,1,ROS,MODET,ROS,CPR120,A104:0,P30:2,200J,360J,360J,50J,75J,90J,ALRT,NZ	48047
99512-000932	CR2-2-000951	LPCR2,WIFI,H,EN-GB,DE,S,1,ROS,ROS,FLX,CPR120,A104:0,P30:2,200J,360J,360J,50J,75J,90J,ALRT,NZ	48047
99512-000933	CR2-2-000952	LPCR2,WIFI,B,EN-GB,DE,A,1,ROS,MODET,ROS,CPR120,A104:0,P30:2,200J,360J,360J,50J,75J,90J,ALRT,NZ	48047
99512-000934	CR2-2-000953	LPCR2,WIFI,H,EN-GB,DE,A,1,ROS,ROS,FLX,CPR120,A104:0,P30:2,200J,360J,360J,50J,75J,90J,ALRT,NZ	48047
99512-000935	CR2-2-000954	LPCR2,CELL,B,EN-GB,DE,S,1,ROS,ROS,FLX,CPR120,A104:0,P30:2,200J,360J,360J,50J,75J,90J,ALRT,NZ	48047
99512-000936	CR2-2-000955	LPCR2,CELL,H,EN-GB,DE,S,1,ROS,ROS,FLX,CPR120,A104:0,P30:2,200J,360J,360J,50J,75J,90J,ALRT,NZ	48047
99512-000937	CR2-2-000956	LPCR2,CELL,B,EN-GB,DE,A,1,ROS,MODET,ROS,CPR120,A104:0,P30:2,200J,360J,360J,50J,75J,90J,ALRT,NZ	48047
99512-000938	CR2-2-000957	LPCR2,CELL,H,EN-GB,DE,A,1,ROS,ROS,FLX,CPR120,A104:0,P30:2,200J,360J,360J,50J,75J,90J,ALRT,NZ	48047
99512-000939	CR2-2-000958	LPCR2,WIFI,B,EN-GB,KO,S,2,ROS,MODET,ROS,CPR120,A104:0,P30:2,200J,360J,360J,50J,75J,90J,ALRT,NZ	48047
99512-000940	CR2-2-000959	LPCR2,WIFI,H,EN-GB,DE,S,1,ROS,MODET,ROS,CPR120,A104:0,P30:2,200J,360J,360J,50J,75J,90J,ALRT,NZ	48047
99512-000941	CR2-2-000960	LPCR2,WIFI,H,EN-GB,DE,A,1,ROS,MODET,ROS,CPR120,A104:0,P30:2,200J,360J,360J,50J,75J,90J,ALRT,NZ	48047
99512-000942	CR2-2-000961	LPCR2,CELL,B,EN-GB,DE,S,1,ROS,MODET,ROS,CPR120,A104:0,P30:2,200J,360J,360J,50J,75J,90J,ALRT,NZ	48047
99512-000943	CR2-2-000962	LPCR2,CELL,H,EN-GB,DE,S,1,ROS,MODET,ROS,CPR120,A104:0,P30:2,200J,360J,360J,50J,75J,90J,ALRT,NZ	48047
99512-000944	CR2-2-000963	LPCR2,CELL,H,EN-GB,DE,A,1,ROS,MODET,ROS,CPR120,A104:0,P30:2,200J,360J,360J,50J,75J,90J,ALRT,NZ	48047
99512-000945	CR2-2-000964	LPCR2,WIFI,H,EN-GB,KO,S,2,ROS,MODET,ROS,CPR120,A104:0,P30:2,200J,360J,360J,50J,75J,90J,ALRT,NZ	48047
99512-000946	CR2-2-000965	LPCR2,WIFI,B,EN-GB,KO,A,2,ROS,MODET,ROS,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,NZ	48047
99512-000947	CR2-2-000966	LPCR2,WIFI,H,EN-GB,KO,A,2,ROS,MODET,ROS,CPR120,A104:0,P30:2,200J,360J,360J,50J,75J,90J,ALRT,NZ	48047
99512-000948	CR2-2-000967	LPCR2,WIFI,B,EN-GB,KO,A,2,ROS,MODET,ROS,CPR120,A104:0,P30:2,200J,360J,360J,50J,75J,90J,ALRT,NZ	48047
99512-000949	CR2-2-000968	LPCR2,CELL,B,EN-GB,KO,S,2,ROS,MODET,ROS,CPR120,A104:0,P30:2,200J,360J,360J,50J,75J,90J,ALRT,NZ	48047
99512-000950	CR2-2-000969	LPCR2,CELL,H,EN-GB,KO,S,2,ROS,MODET,ROS,CPR120,A104:0,P30:2,200J,360J,360J,50J,75J,90J,ALRT,NZ	48047
99512-000951	CR2-2-000970	LPCR2,CELL,B,EN-GB,KO,A,2,ROS,MODET,ROS,CPR120,A104:0,P30:2,200J,360J,360J,50J,75J,90J,ALRT,NZ	48047
99512-000952	CR2-2-000971	LPCR2,CELL,H,EN-GB,KO,A,2,ROS,MODET,ROS,CPR120,A104:0,P30:2,200J,360J,360J,50J,75J,90J,ALRT,NZ	48047
99512-000953	CR2-2-000972	LPCR2,WIFI,B,DE,FR,S,2,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,LU	48047
99512-000954	CR2-2-000973	LPCR2,WIFI,H,DE,FR,S,2,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,LU	48047
99512-000955	CR2-2-000974	LPCR2,WIFI,B,DE,FR,A,2,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,LU	48047
99512-000956	CR2-2-000975	LPCR2,WIFI,H,DE,FR,A,2,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,LU	48047
99512-000957	CR2-2-000976	LPCR2,CELL,B,DE,FR,S,2,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,LU	48047
99512-000958	CR2-2-000977	LPCR2,CELL,H,DE,FR,S,2,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,LU	48047
99512-000959	CR2-2-000978	LPCR2,CELL,B,DE,FR,A,2,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,LU	48047
99512-000960	CR2-2-000979	LPCR2,CELL,H,DE,FR,A,2,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,LU	48047
99512-000961	CR2-2-000980	LPCR2,CELL,B,NL,DE,S,1,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,360J,360J,50J,90J,90J,ALRT,NL	48047
99512-000962	CR2-2-000981	LPCR2,WIFI,B,FR,EN-GB,S,2,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,LU	48047
99512-000963	CR2-2-000982	LPCR2,WIFI,H,FR,EN-GB,S,2,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,LU	48047
99512-000964	CR2-2-000983	LPCR2,WIFI,B,FR,EN-GB,A,2,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,LU	48047
99512-000965	CR2-2-000984	LPCR2,WIFI,H,FR,EN-GB,A,2,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,LU	48047
99512-000966	CR2-2-000985	LPCR2,CELL,B,FR,EN-GB,S,2,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,LU	48047
99512-000967	CR2-2-000986	LPCR2,CELL,H,FR,EN-GB,S,2,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,LU	48047
99512-000968	CR2-2-000987	LPCR2,CELL,B,FR,EN-GB,A,2,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,LU	48047
99512-000969	CR2-2-000988	LPCR2,CELL,H,FR,EN-GB,A,2,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,LU	48047
99512-000970	CR2-2-000989	LPCR2,CELL,B,EN-US,DE,S,1,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,US	48047
99512-000971	CR2-2-000990	LPCR2,WIFI,B,EN-GB,DE,S,1,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,HK	48047

Ref #	Part Number	Description	GMDN
99512-001192	CR2-2-001211	LPCR2,CELL,H,AR,EN-GB,S,2,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,EG	48047
99512-001193	CR2-2-001212	LPCR2,CELL,B,AR,EN-GB,A,2,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,EG	48047
99512-001194	CR2-2-001213	LPCR2,CELL,H,AR,EN-GB,A,2,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,EG	48047
99512-001196	CR2-2-001214	LPCR2,WIFI,B,JA,DE,S,1,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,JP	48047
99512-001197	CR2-2-001215	LPCR2,WIFI,B,JA,DE,S,1,MODET,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,JP	48047
99512-001198	CR2-2-001216	LPCR2,WIFI,B,JA,EN-GB,S,2,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,JP	48047
99512-001199	CR2-2-001217	LPCR2,WIFI,B,JA,EN-GB,S,2,MODET,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,JP	48047
99512-001205	CR2-2-001222	LPCR2,USB,H,EN-GB,DE,A,1,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,AU	48047
99512-001206	CR2-2-001223	LPCR2,USB,B,EN-GB,ZH-TW,S,2,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,AU	48047
99512-001207	CR2-2-001224	LPCR2,USB,H,EN-GB,ZH-TW,S,2,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,AU	48047
99512-001208	CR2-2-001225	LPCR2,USB,B,EN-GB,ZH-TW,A,2,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,AU	48047
99512-001209	CR2-2-001226	LPCR2,USB,H,EN-GB,ZH-TW,A,2,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,AU	48047
99512-001210	CR2-2-001227	LPCR2,USB,B,EN-GB,DE,S,1,MODET,CPR120,A104:0,P30:2,200J,360J,360J,50J,75J,90J,ALRT,NZ	48047
99512-001211	CR2-2-001228	LPCR2,USB,H,EN-GB,DE,S,1,MODET,CPR120,A104:0,P30:2,200J,360J,360J,50J,75J,90J,ALRT,NZ	48047
99512-001212	CR2-2-001229	LPCR2,USB,B,EN-GB,DE,A,1,MODET,CPR120,A104:0,P30:2,200J,360J,360J,50J,75J,90J,ALRT,NZ	48047
99512-001213	CR2-2-001230	LPCR2,USB,H,EN-GB,DE,A,1,MODET,CPR120,A104:0,P30:2,200J,360J,360J,50J,75J,90J,ALRT,NZ	48047
99512-001214	CR2-2-001231	LPCR2,USB,H,EN-GB,KO,S,2,MODET,CPR120,A104:0,P30:2,200J,360J,360J,50J,75J,90J,ALRT,NZ	48047
99512-001215	CR2-2-001232	LPCR2,USB,B,EN-GB,KO,S,2,MODET,CPR120,A104:0,P30:2,200J,360J,360J,50J,75J,90J,ALRT,NZ	48047
99512-001216	CR2-2-001233	LPCR2,USB,B,EN-GB,KO,A,2,MODET,CPR120,A104:0,P30:2,200J,360J,360J,50J,75J,90J,ALRT,NZ	48047
99512-001217	CR2-2-001234	LPCR2,USB,H,EN-GB,KO,A,2,MODET,CPR120,A104:0,P30:2,200J,360J,360J,50J,75J,90J,ALRT,NZ	48047
99512-001218	CR2-2-001235	LPCR2,USB,B,EN-GB,DE,S,1,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,IE	48047
99512-001219	CR2-2-001236	LPCR2,USB,H,EN-GB,DE,S,1,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,IE	48047
99512-001220	CR2-2-001237	LPCR2,USB,B,EN-GB,DE,A,1,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,IE	48047
99512-001221	CR2-2-001238	LPCR2,USB,H,EN-GB,DE,A,1,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,IE	48047
99512-001222	CR2-2-001239	LPCR2,USB,B,EN-GB,DE,S,1,MODET,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,IE	48047
99512-001223	CR2-2-001240	LPCR2,USB,H,EN-GB,DE,S,1,MODET,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,IE	48047
99512-001224	CR2-2-001241	LPCR2,USB,B,EN-GB,DE,A,1,MODET,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,IE	48047
99512-001225	CR2-2-001242	LPCR2,USB,H,EN-GB,DE,A,1,MODET,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,IE	48047
99512-001226	CR2-2-001243	LPCR2,USB,B,EN-GB,DE,S,1,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,GB	48047
99512-001227	CR2-2-001244	LPCR2,USB,H,EN-GB,DE,S,1,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,GB	48047
99512-001228	CR2-2-001245	LPCR2,USB,B,EN-GB,DE,A,1,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,GB	48047
99512-001229	CR2-2-001246	LPCR2,USB,H,EN-GB,DE,A,1,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,GB	48047
99512-001230	CR2-2-001247	LPCR2,USB,B,EN-GB,DE,S,1,MODET,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,GB	48047
99512-001231	CR2-2-001248	LPCR2,USB,H,EN-GB,DE,S,1,MODET,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,GB	48047
99512-001232	CR2-2-001249	LPCR2,USB,B,EN-GB,DE,A,1,MODET,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,GB	48047
99512-001233	CR2-2-001250	LPCR2,USB,H,EN-GB,DE,A,1,MODET,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,GB	48047
99512-001234	CR2-2-001251	LPCR2,USB,B,FR,DE,S,1,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,FR	48047
99512-001235	CR2-2-001252	LPCR2,USB,H,FR,DE,S,1,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,FR	48047
99512-001236	CR2-2-001253	LPCR2,USB,B,FR,DE,A,1,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,FR	48047
99512-001237	CR2-2-001254	LPCR2,USB,H,FR,DE,A,1,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,FR	48047
99512-001238	CR2-2-001255	LPCR2,USB,B,FR,EN-GB,S,2,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,FR	48047
99512-001239	CR2-2-001256	LPCR2,USB,H,FR,EN-GB,S,2,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,FR	48047
99512-001240	CR2-2-001257	LPCR2,USB,B,FR,EN-GB,A,2,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,FR	48047
99512-001241	CR2-2-001258	LPCR2,USB,H,FR,EN-GB,A,2,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,FR	48047

Ref #	Part Number	Description	GMDN
99512-001242	CR2-2-001259	LPCR2,USB,B,ES-ES,DE,S,1,MODET,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,ES	48047
99512-001243	CR2-2-001260	LPCR2,USB,H,ES-ES,DE,S,1,MODET,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,ES	48047
99512-001244	CR2-2-001261	LPCR2,USB,B,ES-ES,DE,A,1,MODET,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,ES	48047
99512-001245	CR2-2-001262	LPCR2,USB,H,ES-ES,DE,A,1,MODET,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,ES	48047
99512-001246	CR2-2-001263	LPCR2,USB,B,ES-ES,EN-GB,S,2,MODET,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,ES	48047
99512-001247	CR2-2-001264	LPCR2,USB,H,ES-ES,EN-GB,S,2,MODET,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,ES	48047
99512-001248	CR2-2-001265	LPCR2,USB,B,ES-ES,EN-GB,A,2,MODET,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,ES	48047
99512-001249	CR2-2-001266	LPCR2,USB,H,ES-ES,EN-GB,A,2,MODET,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,ES	48047
99512-001250	CR2-2-001267	LPCR2,USB,B,EN-GB,DE,S,1,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,AU	48047
99512-001251	CR2-2-001268	LPCR2,USB,H,EN-GB,DE,S,1,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,AU	48047
99512-001252	CR2-2-001269	LPCR2,USB,B,EN-GB,DE,A,1,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,AU	48047
99512-001261	CR2-2-001278	LPCR2,WIFI,B,EN-US,DE,S,1,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,US	48047
99512-001262	CR2-2-001279	LPCR2,WIFI,H,EN-US,DE,S,1,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,US	48047
99512-001263	CR2-2-001280	LPCR2,WIFI,B,EN-US,DE,A,1,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,US	48047
99512-001264	CR2-2-001281	LPCR2,WIFI,H,EN-US,DE,A,1,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,US	48047
99512-001265	CR2-2-001282	LPCR2,WIFI,B,EN-US,ES-ES,S,2,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,US	48047
99512-001266	CR2-2-001283	LPCR2,WIFI,H,EN-US,ES-ES,S,2,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,US	48047
99512-001267	CR2-2-001284	LPCR2,WIFI,B,EN-US,ES-ES,A,2,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,US	48047
99512-001268	CR2-2-001285	LPCR2,WIFI,H,EN-US,ES-ES,A,2,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,US	48047
99512-001269	CR2-2-001286	LPCR2,USB,H,NO,DE,S,1,MODET,FLX,CKB,CPR180,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,NO	48047
99512-001270	CR2-2-001287	LPCR2,USB,H,NO,DE,A,1,MODET,FLX,CKB,CPR180,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,NO	48047
99512-001271	CR2-2-001288	LPCR2,USB,H,FI,DE,S,1,MODET,FLX,CKB,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,FI	48047
99512-001272	CR2-2-001289	LPCR2,USB,H,FI,DE,A,1,MODET,FLX,CKB,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,FI	48047
99512-001273	CR2-2-001290	LPCR2,USB,H,NL,DE,S,1,MODET,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,NL	48047
99512-001274	CR2-2-001291	LPCR2,USB,H,NL,DE,A,1,MODET,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,NL	48047
99512-001275	CR2-2-001292	LPCR2,USB,H,DA,DE,S,1,MODET,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,DK	48047
99512-001276	CR2-2-001293	LPCR2,USB,H,DA,DE,A,1,MODET,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,DK	48047
99512-001277	CR2-2-001294	LPCR2,USB,H,NL,DE,S,1,MODET,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,BE	48047
99512-001278	CR2-2-001295	LPCR2,USB,H,NL,DE,A,1,MODET,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,BE	48047
99512-001279	CR2-2-001296	LPCR2,USB,H,SV,DE,S,1,MODET,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,SE	48047
99512-001280	CR2-2-001297	LPCR2,USB,H,SV,DE,A,1,MODET,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,SE	48047
99512-001281	CR2-2-001298	LPCR2,USB,H,IT,DE,S,1,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,CH	48047
99512-001282	CR2-2-001299	LPCR2,USB,H,IT,DE,A,1,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,CH	48047
99512-001283	CR2-2-001300	LPCR2,USB,H,FR,DE,S,1,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,CH	48047
99512-001284	CR2-2-001301	LPCR2,USB,H,FR,DE,A,1,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,CH	48047
99512-001285	CR2-2-001302	LPCR2,USB,H,DE,EN-GB,S,1,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,CH	48047
99512-001286	CR2-2-001303	LPCR2,USB,H,DE,EN-GB,A,1,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,CH	48047
99512-001287	CR2-2-001304	LPCR2,USB,H,DE,EN-GB,S,1,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,DE	48047
99512-001288	CR2-2-001305	LPCR2,USB,H,DE,EN-GB,A,1,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,DE	48047
99512-001289	CR2-2-001306	LPCR2,USB,H,DE,EN-GB,S,1,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,LI	48047
99512-001290	CR2-2-001307	LPCR2,USB,H,DE,EN-GB,A,1,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,LI	48047
99512-001291	CR2-2-001308	LPCR2,USB,H,DE,EN-GB,S,1,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,AT	48047
99512-001292	CR2-2-001309	LPCR2,USB,H,DE,EN-GB,A,1,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,AT	48047
99512-001293	CR2-2-001310	LPCR2,USB,H,IT,DE,S,1,MODET,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,IT	48047

Ref #	Part Number	Description	GMDN
99512-001382	CR2-2-001399	LPCR2,WIFI,B,EN-GB,AR,A,2,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,JO	48047
99512-001383	CR2-2-001400	LPCR2,WIFI,H,EN-GB,AR,A,2,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,JO	48047
99512-001384	CR2-2-001401	LPCR2,CELL,B,EN-GB,DE,S,1,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,JO	48047
99512-001385	CR2-2-001402	LPCR2,CELL,H,EN-GB,DE,S,1,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,JO	48047
99512-001386	CR2-2-001403	LPCR2,CELL,B,EN-GB,DE,A,1,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,JO	48047
99512-001387	CR2-2-001404	LPCR2,CELL,H,EN-GB,DE,A,1,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,JO	48047
99512-001388	CR2-2-001405	LPCR2,CELL,B,EN-GB,AR,S,2,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,JO	48047
99512-001389	CR2-2-001406	LPCR2,CELL,H,EN-GB,AR,S,2,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,JO	48047
99512-001390	CR2-2-001407	LPCR2,CELL,B,EN-GB,AR,A,2,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,JO	48047
99512-001391	CR2-2-001408	LPCR2,CELL,H,EN-GB,AR,A,2,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,JO	48047
99512-001392	CR2-2-001409	LPCR2,WIFI,B,EN-GB,DE,S,1,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,BH	48047
99512-001393	CR2-2-001410	LPCR2,WIFI,H,EN-GB,DE,S,1,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,BH	48047
99512-001394	CR2-2-001411	LPCR2,WIFI,B,EN-GB,DE,S,1,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,BH	48047
99512-001395	CR2-2-001412	LPCR2,WIFI,H,EN-GB,DE,S,1,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,BH	48047
99512-001396	CR2-2-001413	LPCR2,WIFI,B,EN-GB,DE,A,1,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,BH	48047
99512-001397	CR2-2-001414	LPCR2,WIFI,H,EN-GB,DE,A,1,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,BH	48047
99512-001398	CR2-2-001415	LPCR2,WIFI,B,EN-GB,AR,S,2,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,BH	48047
99512-001399	CR2-2-001416	LPCR2,WIFI,H,EN-GB,AR,S,2,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,BH	48047
99512-001400	CR2-2-001417	LPCR2,WIFI,B,EN-GB,AR,A,2,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,BH	48047
99512-001401	CR2-2-001418	LPCR2,WIFI,H,EN-GB,AR,A,2,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,BH	48047
99512-001402	CR2-2-001419	LPCR2,CELL,B,EN-GB,DE,S,1,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,BH	48047
99512-001403	CR2-2-001420	LPCR2,CELL,H,EN-GB,DE,S,1,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,BH	48047
99512-001404	CR2-2-001421	LPCR2,CELL,B,EN-GB,DE,A,1,ROS,MODET,ROS,FLX,CPR120,A104:0,P30:2,200J,300J,360J,50J,75J,90J,ALRT,BH	48047
99512-001405	CR2-2-001422	LPCR2,CELL,H,EN-GB,DE,A,1,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,BH	48047
99512-001406	CR2-2-001423	LPCR2,CELL,B,EN-GB,AR,S,2,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,BH	48047
99512-001407	CR2-2-001424	LPCR2,CELL,H,EN-GB,AR,S,2,ROS,MODET,ROS,FLX,CPR120,A30:2,P30:2,200J,300J,360J,50J,75J,90J,ALRT,BH	48047



Signed April 15, 2019



Caleb Simpson
Director | Regulatory and Clinical Affairs

This declaration is issued under the sole responsibility of the manufacturer.

This declaration applies to CE marked devices produced after the date of issuance of this declaration and before it is either superseded by another declaration or withdrawn.