

Filtering Half Mask

How to use



1. Check the mask is worn
2. Place the mask over the nose and mouth
3. Pull the top strap over the head
4. Adjust the bottom strap

Notes

1. Check the mask is worn
2. Place the mask over the nose and mouth
3. Pull the top strap over the head
4. Adjust the bottom strap



Produced in
Germany

Filtering Half Mask

EN149:2001+A1:2009 FFP2 NR



6 973149 410037

Soft and Breathable

Safe and Reliable

CE 2834



Waterproof
Temperature
Instructions

EN149:2001+A1:2009 FFP2 NR
Produced in Germany
with a maximum efficiency of 94%
and a maximum flow rate of 120 l/min
at a pressure of 250 Pa. The shelf life is 2 years from the date of production.
Check the instructions for use.
Detailed information: EN149:2001+A1:2009
Production date: see the instructions

30 PCS