

FITAGE

BLOOD PRESSURE MONITOR



2 Users



Large LCD Display



Audio Function



ESH Approved



Heart Rate Monitor



300 Memory Slots

EASY TO USE

STEP 1

Wrap cuff around upper arm



STEP 2

Press "START/STOP" button to start measurement



STEP 3

Memory check
Press "M" button to view previously recorded measurements



PRODUCT SPECIFICATIONS:

Product Description	FITAGE Arm-type Blood Pressure Monitor	
Model	FITAGE-06	
Display	LCD Digital Display	
Measurement Method	Oscillometric Method	
Pressurization	Automatic Pressurization	
Memory	120 Memories in Two Groups with Date and Time	
Power Source	4x AA batteries and USB-C cable	
Operating Environment	Temperature	10°C ~ 40°C (50°F~104°F)
	Humidity	15% ~ 93%RH
	Pressure	700hPa ~ 1060hPa
Storage Environment	Temperature	-25°C~70°C (-13°F~158°F)
	Humidity	≤93% RH
Classification	Internal Powered Equipment, Type BF, Cuff is the Applied Part	
Ingress Protection Rating	IP20, Indoor Use Only	

FITAGE B.V.,
Kabelweg 80, 2371 DX
Roelofarendsveen, Nederland
support@fitage.nl



Shanghai International Holding Corp. GmbH (Europe)
Eiffestrasse 80, 20537 Hamburg, Germany



Designed by FITAGE in Europe, Assembled in China

FITAGE

BLOOD PRESSURE MONITOR

