



Your praise is our reward!
Like our Oximeter? Spread the vibes.
Leave your review.

We depend on less than 2% of our customers
that give 1 minute of their
time to leave a review. Please be one
of them! Thank you so much!

Instruction
video link:
lead.me/OXI-27



SCAN ME

or
scan the QR-code

UK REP: SUNGO Certification Company Limited
Unit G25 Waterfront Studios 1 Dock Road,
London United Kingdom, E16 1AH

MedPain GmbH
Mies-van-der-Rohe-Strasse 8,
8087 Munich, Germany

Shenzhen Juniper Medical Equipment Co., Ltd.
Address: D Building, No. 71, Xixian Road, Fuyong Street,
Baoan, Shenzhen, Guangdong, China

Made in China

CE 0482

Manufactured for Proven
Burg de Raadstingel 71
3311 JG Dordrecht
packaging V1.1

8 719699 636385

For Sports &
Aviation use only



Battery
%SpO2 or
oxygen saturation
Pulse rate
intensity graph
Pulse rate
per minute

It's our passion to bring professional quality right at your doorstep.
So you can monitor your health with devices you can trust.
We want to learn, grow and improve. That's why we welcome your feedback.
We're listening.

www.proven.com

Keep dry.
Upward

Non-condensing
RH:10%-92%

Temperature
T:0-42°C

Atmospheric pressure
50 kPa-106 kPa

Our Passion

Fingertip Pulse Oximeter

- ✓ Bar-Graph for Pulse Intensity Model: OXI-27iWB
- ✓ Automatic Pulse Detection
- ✓ One Click Measurement



One Click
Measurement



For SpO₂
Saturation



Bar-Graph for
Pulse Intensity