

## Product parameters

Product Name : V Face Beauty Meter

Net product : 200g

Product Model : HY18

Charging mode : USB

Product Size : 165×118

Rated voltage : 3.7V

Massage time : 10 minutes

Packing list: Host lift belt, USB charging line, instruction manual, packing box

Maintenance and cleaning : if there is dirt in the appearance of the product, clean with dry, unscrewed cloth or paper towel, no need to use detergent, do not wash!

Note : do not overuse, if there are adverse reactions, it is recommended to stop using. When the electricity is insufficient, charge in time to avoid affecting the life of the product. It is forbidden to disassemble or repair the product. placed where children can not touch.