



EASY ICE

READY-TO-USE, INSTANT ICE
INDICATED FOR HAEMATOMAS,
SPRAINS AND INJURIES

Dispotech disposable instant ice bag, available in 2 versions, made in polyethylene (PE), or made in non woven (TNT) used for sensible part of body; is suitable in case of haematomas, sprains, injuries, hyperthermia and whenever cold applications is required.

EASY TO USE

Break the inside water bubble, by squeezing the bag with both hands. Then shake it in order to mix the contents. Apply on the interested part.

Medical device class IIa



Dispotech Srl

Via del piano 29, 23020 Gordona (SO) Italy
Tel +39 0343/36711, Fax +39 0343/36567
e-mail: info@dispotech.it - www.dispotech.com



dispotech.com