

REF 14929

**Knee Strap**  
**Knee Strap**

**Hansaplast**  
**SPORT**

**XX** This Hansaplast Knee Strap is designed to help relieve symptoms of Patello-femoral Pain Syndrome and tendonitis of the knee.

**Helps provide relief from kneecap pain**

The built-in pressure strip provides targeted support to the tendons below the kneecap to guide the kneecap and helps relieve the associated knee pain

Easily adjustable strap for secure custom fit

Soft-touch, latex-free, moisture releasing material for extra comfort

**For left and right knee, washable and reusable**  
**For left and right knee, washable and reusable**

Wear during activities that lead to discomfort

- Position the pressure pad directly under kneecap, facing skin.
- Loop strap through buckle.
- Adhere fastener to required fit.

**Seek medical advice if in any doubt about the injury.**  
Store in a cool, dry place.

**Hansaplast**  
**SPORT**

**Knee Strap**  
**Knee Strap**

**Helps provide relief from kneecap pain**  
**Helps provide relief from kneecap pain**

**Beiersdorf**  
Beiersdorf AG  
20245 Hamburg, Germany  
Hansaplast / Hansaplast Logo  
= reg. tm. of Beiersdorf AG, Germany

**Made in Taiwan**

Materials: Polyester, Nylon, Ethylene Vinyl Acetate, Nitrile Rubber, Polyethylene

**www.hansaplast.com**

CE

**1 Bandage**      **1 Bandage**

EAN 13 SC2 position only (100% magnification)  
Minimum dimensions:  
Bar width 31.35mm  
Short bar height 22.85mm  
Left light margin 3.63mm  
Right light margin 2.31mm

Exp. Dat. / LOT

DRVxxxx.xx

**Hansaplast**  
**SPORT**

For more than 100 years Beiersdorf offers cutting-edge plaster innovations and Hansaplast is the No. 1 choice for millions of people all around the world. Our experience and expertise in plaster and tapes technology combined with Beiersdorf's extensive skin care competencies enables us to develop products of highest quality and effectiveness to keep you active and moving.

For more than 100 years Beiersdorf offers cutting-edge plaster innovations and Hansaplast is the No. 1 choice for millions of people all around the world. Our experience and expertise in plaster and tapes technology combined with Beiersdorf's extensive skin care competencies enables us to develop products of highest quality and effectiveness to keep you active and moving.

**Adjustable / Adjustable**  
31,8 - 44,5cm

**SUPPORT LEVEL - SUPPORT LEVEL**  
Moderate / Moderate  
**FIRM / FIRM**  
Extra Firm / Extra Firm

**1 Bandage**

| MASTERARTWORK           |                                  |  |
|-------------------------|----------------------------------|--|
| BRV                     | 10123704.01                      |  |
| Product Name            | Hansaplast Sport Knee Strap      |  |
| Packaging               | Foldingbox                       |  |
| Agency                  | DSN GmbH                         |  |
| Date                    | 24.06.14 MHA_DSN No.: 11BED00131 |  |
| Changes to this Artwork | New Masterartwork                |  |
| Optional text           | Storage conditions               | Attention: Declare in Middle East countries (e.g. UAE, Saudi Arabia, Oman, Kuwait, Qatar, Lebanon) "Store in a dry place not above 25°C".  |
|                         | Distributed by                   | In case distribution address(es) will be declared then add "distributed by" before it.   |
|                         | Legal Manufacturer               | Always declare this symbol in front of the Legal Manufacturer address - except in cases not sufficient printing area is available:<br>Beiersdorf AG<br>20245 Hamburg, Germany                                      |
|                         | Made in Taiwan                   | Declaration "Made in Taiwan" is obligatory for countries outside of the EU. For countries being part of the EU the declaration "Made in Taiwan" is not obligatory - in this case please check locally if required. |

Job Date: 19.05.14 MHA  
 1. correction: 20.05.14 ZK  
 2. correction: 27.05.14 ZK  
 3. correction: 27.05.14 MHA  
 4. correction: 11.06.14 MHA  
 5. correction: 13.06.14 MHA  
 6. correction: 19.06.14 MHA  
 7. correction: 20.06.14 MHA  
 8. correction: 24.06.14 MHA

| Colors / decoration | Printing technology/ finishing technology |
|---------------------|---|
| 1. PANTONE 2747 C   | Offset                                    |
| 2. PANTONE 186 C    | Offset                                    |
| 3. Cyan             | Offset                                    |
| 4. Magenta          | Offset                                    |
| 5. Yellow           | Offset                                    |
| 6. Black            | Offset                                    |

**Attention**  
Die cut for control only and not for production. Please use current die cut.

|                 |  |   |
|-----------------|--|---|
| Additional text |  | The green dot symbol is a third party's trademark. Please check outside of Germany if allowed and use it only in those countries where you have concluded a respective license agreement with the owner of the mark. This has not been done centrally. Otherwise please delete this symbol and substitute it by your locally required recycling symbol, if any. |
|                 |  | Please check locally if these two recycling symbols are allowed and/or required in your country.  |
|                 |  | In local artworks the certification number of the current printer needs to be checked and included.   |

95,25 mm      38,1 mm      95,25 mm      37,3 mm