





circulation
maxx
reviver

bioenergiser



circulation
maxx
reviver


Programme timer: 1 to 60 mins



14 Electrical Muscle Stimulation
pattern cycles for foot plates



10 Electrical Muscle Stimulation
pattern cycles for TENS pads



Infrared, hand-held remote control



2 pairs of TENS electro pads



EC REP

DISTRIBUTED BY:
Orange Planet
De Els 3A, 5107 RK Dongen
The Netherlands

PRODUCT OWNER:
Bioenergiser Ltd.
Unit 2 Frederick Road, Hoo Farm Estate
Kidderminster, Worcestershire DY11 7RA

Contraindications: DO NOT USE if you have a Pacemaker, current DVT or cardiovascular problems, in pregnancy or with epilepsy. 18 years old + only. Not suitable for children.



© Bioenergiser™ Circulation Reviver™ are trademarks of
Bioenergiser (UK) Ltd. All rights reserved.
Manufactured for Bioenergiser (UK) Ltd. in the PPR

ANOTHER GREAT PRODUCT DISTRIBUTED BY

