


Designed by A&D Medical, Tokyo, Japan.  
 Made in China  
  
  
 0434  
 A&D Company Limited  
 3-23-14 Higashi-Kobukuro, Toshima-ku Tokyo 170-0013 JAPAN  
<http://www.aandd.jp>

 Caution Please read instruction manual thoroughly before using.


- Soft, flexible tip for safer measuring
- Mercury-free
- 60 second reading
- Measurement time will vary by individual and method used.
- For Oral, Underarm or Rectal use
- Beeps when measurement complete
- Last temperature reading memory
- Accurate readings within  $\pm 0.1^{\circ}\text{C}$
- Water-resistant sensor
- Auto power off approximately 10 minutes
- $^{\circ}\text{C}/^{\circ}\text{F}$  conversion

**Features**

# Digital Thermometer



## Digital Thermometer



- 60 second reading
- Accurate readings within  $\pm 0.1^{\circ}\text{C}$

**With Flexible tip**

  
 Instruction manual in 5 languages

