

Android download



①:Open the **Google Play**



②:Search **EndscopeTool** and click to download



③:click to download **EndscopeTool**




④:Click to open **EndscopeTool**

Apple download



①:Open the **App Store**



②:Search **EndscopeTool** and click to GET



③:Open **EndscopeTool**



2 9 5 0 0 2 1 0 9 3 5 4 4
URGOODS Oorsmee...reiniger LVB

Product Name: Earwax Removal Tool

Network Standard: IEE802.11b/g/n

Working Temperature: -10~40°C

Executive Standard: GB4943

Product Model: Y8

Battery Capacity: 130mAh

Working Time: 30 minutes

Input Current:DC 5V300m A



Made in China