

MCS FORMULAS EGCG-80% GREEN TEA,

contains high dose [80%] EGCG [Epigallocatechin gallate], one of the most abundant catechin in green tea. EGCG is a polyphenol intensively researched for its potential to affect human health and disease. To maximize its supportive role, this product has been designed to contain one of the largest amounts of EGCG available in one capsule [320 mg], decaffeinated. The high concentration allowed us to maintain strength and reduce capsule size.

If you are pregnant, nursing, taking any medications, please consult your healthcare practitioner before taking any dietary supplements. Do not use if inner bottle seal is missing, torn or damaged in any way. This product is not intended to diagnose, treat, cure, or prevent any disease.

Food supplements should not be used as a substitute for a varied and healthy lifestyle. Keep out of reach from children. Do not exceed the daily intake recommended by your healthcare provider.

Produced with high quality standards, and packaged in Europe.

Carefully manufactured for:
MCS Formulas B.V.
Postbus 78, 5600 AB
Eindhoven, The Netherlands

Contact us for more information:
www.mcsformulas.com



49.5 g e



We donate 50% to improve
and extend life of cancer patients
Learn more on www.mcsformulas.com



— HIGH STRENGTH —

EGCG-80%

400 mg Green Tea Extract in One Capsule

High Dose & High Purity
Decaffeinated

Dietary Supplement

100
VEG. CAPSULES

NON
GMO

GLUTEN
FREE

SOY
FREE

SUPPLEMENT FACTS

SERVING SIZE: 1 CAPSULE
SERVINGS PER CONTAINER: 100

Amount Per Capsule		%DV
Green Tea Extract [Camellia sinensis], Leaf	400 mg	*
Std. for a minimum of 80% Epigallo-catechin-3 gallate [EGCG]	320 mg	*

*Daily Value [DV] not established

OTHER INGREDIENTS: HYPROMELLOSE [VEGAN CAPSULE], MICROCRYSTALLINE CELLULOSE.

DOES NOT CONTAIN: GMO, GLUTEN, WHEAT, DAIRY, YEAST, SUGAR, SODIUM, ARTIFICIAL FLAVOR, SWEETENER, PRESERVATIVES & COLOR. NOT IRRADIATED. SUITABLE FOR VEGANS.

SUGGESTED USE: AS A DIETARY SUPPLEMENT, ADULTS TAKE 1 CAPSULE DAILY 30 MIN. BEFORE FOOD OR AS DIRECTED BY A HEALTHCARE PROFESSIONAL.

STORE: IN A COOL DRY PLACE.

BEST BEFORE: SEE BASE
BATCH NUMBER: SEE BASE
PRODUCT CODE: EGC-2-MCS

